


6' WORKBENCH TWS-1900

Teclab workbenches share common design features that have earned for Teclab the reputation in the industry as the "professional's bench." Typical features include an abundance of grounded outlets within easy reach; an infinite variety of drawer and door combinations to meet every job and storage requirement; expansive work surfaces, shelves, and back guards in a choice of task-oriented materials.



FEATURES

One **C-1521-EP-1** three drawer, 15" wide Base Cabinet with 5 grounded outlets, plus illuminated combination rocker switch

One **C-3010-EP-1** two drawer, two door, 30" wide Base Cabinet with 10 grounded outlets, plus illuminated combination rocker switch

One **FR6-25** 25" Foot Rest

One **APD-2200** 22" wide Apron Drawer

Two **AF-01.5** 1.5" Apron Fillers

Two **SS-16-EP-1** 16" tall Shelf Supports with 5 grounded 15A outlets, plus illuminated combination rocker switch.


One **TP1-3072-PVC** 30" x 72" x 1 1/4" plastic laminate Work Surface

One **SP1-1272-PVC** 12" x 72" x 1 1/4" Shelf of same material

Two **BGP-72-PVC** 4" x 72" x 1 1/4" Back Guards of same material

FINISH OPTIONS


Powder coated steel components.



Custom paint color matching is available. Stainless Steel is also available.

WORK SURFACE COLOR OPTIONS

Visual comparison. Color choices vary depending on work surface type.



PD - TWS-1900 - v.1.10