



LABORATORY AND TECHNICAL FURNITURE



6450 Valley Industrial Drive
Kalamazoo, MI 49009

269-372-6000
info@teclab.com
teclab.com

TECLAB IS RECOGNIZED AS A LEADER IN THE DESIGN, MANUFACTURE AND INSTALLATION OF LABORATORY FURNITURE, ERGONOMIC WORKSTATIONS, WORKBENCH SYSTEMS, AND STORAGE EQUIPMENT.

Whether you're looking to create a single workstation or an entire laboratory facility ,Teclab can supply all the essential components for a highly efficient and productive work space.

We are an American company that's been in the laboratory furniture business for more than 40 years supplying top companies in the energy, pharmaceutical, healthcare, biotech and research fields as well as government and educational facilities. These organizations rely on Teclab to supply the best products and workmanship.

Customer service has always been a top priority at Teclab and our extremely knowledgeable staff can also assist you in planning your laboratory. From room layout to installation, our in-house designers and engineers use their decades of laboratory design experience to customize your application needs.

You can call us at **269.372.6000** and speak one on one with a helpful and friendly member of Teclab's planning and design team or send us an email at **info@teclab.com** and start creating your new laboratory.

Office Hours:

8:00 a.m. - 5:00 p.m. EST

We're here to answer any questions you may have. It's our goal to make your experience with Teclab a pleasant and rewarding one that will also make your laboratory the best it can be.


FREE PLANNING AND DESIGN

LET US DO IT FOR YOU!

FREE!

Call us today at

269.372.6000



IN GOOD COMPANY



TESLA

American Airlines



SPACEX

Johnson & Johnson



NASA

BorgWarner

BLUE ORIGIN

COLGATE-PALMOLIVE



ExxonMobil

Los Alamos
NATIONAL LABORATORY



BOEING

GENERAL DYNAMICS



3M

LOCKHEED MARTIN

NOKIA

Bell Laboratories



KIA

BOSCH



FedEx

NORTHROP GRUMMAN

xerox

AIRBUS



SONY

BASF

HONDA

TOYOTA

IBM

HYUNDAI

LABORATORY CASEWORK COMPONENTS, ACCESSORIES AND FUME HOODS

At Teclab, today's lab comes with the built-in adaptability you'll need to change tomorrow.

Highest quality. Best built. Most durable. Most versatile. Take your pick. They're all descriptions that have become synonymous with Teclab in the last 30 years. But our fanatical commitment to building the best is more than just talk.

The extraordinary fit and finish of every Teclab product reflects that devotion to providing what customers want and need. Customers see our quality in clean, functional lines. They hear it in the sound of a drawer sliding effortlessly closed. The feel it in the solidity of a well made unit. And they experience it in the attention to the smallest of details of design assistance, installation and after-sale service.


Teclab builds modular cabinetry of superior strength, thanks to twin full-perimeter frames. We also set door and drawer faces into the cabinet shell to protect the contents from spills. Specify 18-, 24-, 30-, 36-, 42- or 48-inch widths.

- Base cabinets, storage cabinets and components in a choice of stainless steel or steel painted in one of 12 standard or any custom-matched color.
- Wide range of work surfaces, depending on lab requirements.
- Accessories, including under-shelf lights, task lights, sinks, chairs, drying racks, service fixtures, bin rails, keyboard trays, monitor arms, CPU holders and tack boards.
- Hands-on lab planning and design assistance from Teclab application engineers.
- Fume hoods and enclosures.



The Base Cabinet is the foundation upon which all Teclab equipped laboratories are created.

Teclab's commitment to quality design can be seen in the details of our cabinets. We start with furniture grade cold rolled steel and use only the finest materials available throughout the manufacturing process. The results are components that meet or exceed the requirements of today's high-tech labs and manufacturing facilities.




SEE OUR FULL SELECTION OF BASE CABINETS AND OPTIONS IN THE LATER PAGES OF THIS CATALOG.

BASE CABINET COMPARISON


Teclab Base Cabinets are offered in two distinct styles.

LABORATORY BASE CABINETS

Teclab's Laboratory(Wetlab) style cabinets are 22" deep. When installed under a standard 30" countertop, there is ample room for plumbing and other utilities behind the cabinets. The drawer and door faces of Lab Cabinets are inset into the cabinet to protect the contents from spills.




BENCH HEIGHT LAB CABINET




DESK HEIGHT LAB CABINET

TECHNICAL BASE CABINETS

Teclab's Technical(TWS) style cabinets are 28" deep. The extra depth allows for massive amounts of storage and increased stability in free-standing installations. TWS Cabinets are often used as a basis for creating a Workbench.



BENCH HEIGHT TWS CABINET



DESK HEIGHT TWS CABINET

Our Hoods incorporate a sleek exterior with a molded one-piece fiberglass liner. The one-piece liner of specially-formulated, fiberglass-reinforced polyester offers corrosion and fire resistance and easy clean up.

Laboratory Hoods are offered in 4', 5', 6' and 8' widths so their dimensions exactly match standard laboratory casework. Choose from hoods designed to be used with a remotely-located blower or hoods that include a built-in blower.

Models with built-in blower are available. Quiet built-in blower is belt-driven with molded thermoplastic housing and coated aluminum impeller that is non-sparking and corrosion-resistant. Available with explosion-proof or standard motor.

LED lighting illuminates the interior. Located outside the hood interior for corrosion-resistance and easy replacement. Explosion-proof models are supplied with an incandescent light fixture.

Color-coded service fixtures for gas, air, water, vacuum and other services have remote controls for use regardless of the sash position. Two service fixtures are pre-plumbed on fixtured models. Each hood is factory prepared for up to eight service fixtures (four on each side).

Molded Epoxy-resin Work-surface. Our standard fume hoods include the required work-surface.*

Performance tested to ASHRAE 110-1995.

Front and side panels may be easily removed for lamp replacement and access to electrical supply connections.

Service access panels allow accessibility to plumbing from the front of the hood.

Ergonomic air foil allows air to sweep the work surface for maximum containment. Airflow openings pull inflow air from under the air foil so that clean air continually flows over the air foil creating a constant barrier of protection from contaminants.

Durable and attractive exterior is glacier white, dry powder epoxy-coated steel.

Pre-set baffle provides a uniform face velocity pattern across the sash opening of the hood.

One-piece molded fiberglass liner gives superior corrosion and chemical resistance, durability and light reflectivity. Its smooth, radiused corners make cleaning easy.

*Custom sized fume hoods and specialized fume hoods will require a custom work-surface tailored to suit the environment.



Laboratory Hoods for use with remote blower

4 Feet, part# FH-48-EPLBSCSTGA

5 Feet, part# FH-60-EPLBSCSTGA

6 Feet, part# FH-72-EPLBSCSTGA

8 Feet, part# FH-96-EPLBSCSTGA

Laboratory Hoods with built-in blower

4 Feet, part# FH-48-EPLBSCSTGABB

5 Feet, part# FH-60-EPLBSCSTGABB

6 Feet, part# FH-72-EPLBSCSTGABB


Molded Epoxy-resin Work-surface. Our standard fume hoods include the required work-surface.*



CANOPY HOODS

Canopy hoods can be built to any dimension needed. Typically, the depth is 24" and the width varies to suit the space.

Canopy Hoods include a heavy-duty wall mounting bracket system for simple, solid installation.




EXTRACTION ARMS

The ideal extractor for laboratory environments.

With its optimal design, the Ø 3 inch Extraction Arm has a very low pressure drop, which provides many valuable benefits:

- Low pressure drop saves energy.
- Air flow noise is reduced.
- Lower pressure drop is achieved without selecting a larger diameter extractor.
- Lower pressure drop allows the extraction arm to be combined with additional extraction systems.



AIRFLOW MONITORS

Airflow monitors sense and alert the operator when the airflow of a fume hood is outside of normal range.



The Industry Standard!

Teclab offers the standard in workbench design and construction. Component for component, Teclab furniture has been designed to offer you the ultimate in durability, flexibility, and performance. Whether it's the 86 critical weld locations or any of the many unique Teclab design features, you'll be impressed by Teclab's commitment to detail. When you demand quality design, specify Teclab.

Teclab's workbench product line consists of several hundred cabinet designs, wall and storage cabinets, height adjustable bench frames, electrical raceways, overhead lighting, work-surfaces, and ergonomic seating. Teclab offers a complete line of Electrostatic Dissipative (ESD) workbenches and accessories to protect sensitive electronic devices.



6' WORKBENCHES

TWS-1000



TWS-1020



TWS-1030



TWS-1100



TWS-1200



TWS-1300



TWS-1400



TWS-1500



TWS-1510

6' WORKBENCHES**TWS-1600****TWS-1700****TWS-1800****TWS-1900****TWS-2000****TWS-3000****TWS-3010****TWS-3012****TWS-3014**

6' WORKBENCHES**TWS-3016****TWS-3018****TWS-3020****TWS-4000****TWS-5000****TWS-6000****TWS-7000****TWS-8000****TWS-9000**

8' WORKBENCHES



TWS-2100



TWS-2200



TWS-2300



TWS-2400



TWS-2500



TWS-2600



TWS-2700



TWS-2800



TWS-4540



TWS-4550



TWS-4560

10' WORKBENCHES**TWS-4600****TWS-4700****TWS-4800****12' WORKBENCHES****TWS-3600****TWS-3700****TWS-3800****TWS-3900****TWS-4100****TWS-4200****TWS-4300****TWS-4400****TWS-4500**

CORNER WORKBENCHES



TCU-1500



TCU-1515



TCU-1520



TCU-1522



TCU-1530



TCU-2222

DESK HEIGHT WORKBENCHES



TWS-5100



TWS-5200



TWS-5300



TWS-5400

ESD OPTIONS



A Teclab Static Control Bench is a primary component of any established ESD Control program.

ESD (Electro-Static Discharge)


The sudden transfer of free electrons from a positively charged object to a negatively charged object.

ESD events occur far more frequently than we realize, and many ESD events go unnoticed. Approximately 90% of all ESD events are generated by the human body, but in order for the human body to feel static electricity, the discharge must reach at least 3,500 volts.


Something as simple as the movement of the technician at the workbench, however, can generate as many as 6,000 volts of static. And costly damage to some electronic devices can occur with as little as 15 volts of static. The unknown nature of ESD translates to an invisible danger in the form of product failures, quality problems and unscheduled downtime.




ESD Wrist Strap



ESD Drawer Mat




Common Point Ground System



ESD Ground Jacks




Isolated Components



ESD Floor Mat




Constant Ground Monitor



ESD Bench Mat



ESD Seating



ESD Work Surface

HEAVY DUTY INDUSTRIAL WORKBENCH

Teclab's TES Line showcases height adjustable and fixed height workbenches with steel uprights supporting overhead lights, storage components and many bench accessories to keep your production components close at hand. We'll custom-make your idea of the perfect bench with the same attention to detail we use in building the world's finest electronic workstations. In the end, you'll receive the superior-quality, high-performance bench you want and need.

CONFIGURATIONS AND OPTIONS

1 Overhead Rack

Accommodates anything from a light to a tool trolley or compressed-air line.

2 Tool Trolley

Keeps assemblers working efficiently. Hook up and go.

3 Adjustable Shelf

Puts tools and parts close at hand, yet keeps them out of the way.

4 Outlet Strip

Puts the power where you need it, in any configuration or amperage.

5 ESD Wrist Strap

Helps eliminate ESD damage caused by human handling

6 Height Adjustment

Maintain the right ergonomics for individual workers; choose from three height adjustment options, including electric, mechanical or manual.

7 Foot Rest

A standard option. The Foot Rest is a structural element and a sensible, simple way to foster greater worker comfort.

8 Steel Overhead Cabinet with Flipper Door

Convenient storage that keeps clutter off the work surface. Available in an array of sizes and colors.

9 LED Light Kit

Under Shelf Light Kits install under Shelves or Cabinets and provide diffused light to the work surface.

10 Uprights

Allow use of multiple accessories, including cabinets, shelves, parts bins, and outlet strips.

11 Parts Bins

The easy way to keep small parts an arm's length away.

12 Work Surface

Choose from a variety of materials to fit your bench and need.

13 Steel Hanging Cabinet

Double-wall construction, a variety of drawer configurations, recessed handles and smooth gliding drawer slides ensure ease of use and maximum durability.

14 Back Support Panel

A standard option. A standard option. The Panel adds strength and stiffness.

White Dry Erase Board or Steel Pegboard

The perfect place to draw notes, or hang tools.

ESD

All Teclab Workbenches can be specified to conform to ESD Standards.


Mobile

Swivel and Rigid Casters with brakes and locks available.


CPU Holders, Keyboard Trays, Monitor Arms.

STANDARD FEATURES

- Durable textured powdercoat finish
- Extraordinary Teclab quality
- Teclab 10-year warranty
- Made in USA
- Adjustable Leveling Feet
- Weight Capacity: 750 lbs.




STANDARD DIMENSIONS



* HEIGHT SHOWN IS BASED ON A 1 1/4\" THICK WORK SURFACE. OVERALL HEIGHT VARIES DEPENDING ON THICKNESS.
* BENCHES LONGER THAN 72\" REQUIRE A CENTER SET OF LEGS.

FINISH OPTIONS


Powder coated steel components.



Custom paint color matching is available. Stainless Steel is also available.

WORK SURFACE COLOR OPTIONS

Visual comparison. Color choices vary depending on work surface type.



VIVI by Teclab.

VIVI Cabinets are constructed of only top quality, heavy-duty, furniture-grade, cold-rolled stainless steel, double-braced and welded at all critical joints; twin full-perimeter frames formed of heavy gauge stainless steel, welded to provide strength and rigidity; heavy-duty drawers that glide effortlessly on precision steel ball bearings; cabinet doors that open and close with authority and precision using self-closing hinges.



Optional Kick-plate snaps into position and keeps debris in check.



Elegant interiors contain no obstructions providing easy to clean hygienic surfaces.



VIVI Cabinets and components lock together in perfect alignment using our Tec-Lock system.

The Teclab Industrial Storage Cabinet. Setting a new standard in storage!

That's exactly what you'll be thinking when you take a good look at the new Teclab Industrial Storage Cabinet. The Teclab design team has redefined the storage cabinet with Teclab's industry leading standard of quality and durability.

The Industrial Storage Cabinet features stronger construction with heavier material than competitive cabinets. With exclusive Teclab features including rugged internal frames, reinforced drawers with smooth operating slide mechanisms and easy to use four-point in cabinet leveling.

**FEATURES**

- 450 lb. drawer slide capacity
- 4800 lb. cabinet capacity (based on caster capacity)
- Full extension drawer access
- Integrated drawer divider system within every drawer
- Fully powder coated components to resist rust and provide a durable scratch-resistant surface
- Teclab 10-year warranty

FRAME

- Fully welded frame
- 12 ga. horizontal frames used in such a way as to reestablish a standard in cabinet rigidity negating the need for a separate base
- 16 ga. vertical frames punched in 50 mm increments allowing extensive drawer sizing options

CABINET OPTIONS

- Available in 18", 24", 30", 36", 42" and 48" widths
- Numerous drawer heights available
- Optional central lock system
- Optional Interlocking one drawer only access
- Available in many colors and stainless steel

LOWER END OPTIONS

- Forklift accessible legs
- Mobile casters
- Removable skirts
- Also available without legs for mounting to a platform or to fit within a tight space



The Teclab Storage cabinet. Setting a new standard in storage!

That's exactly what you'll be thinking when you take a good look at the new Teclab storage cabinet. The Teclab design team has redefined the storage cabinet with Teclab's industry leading standard of quality and durability.

The Storage cabinet features beefier construction with heavier material than competitive cabinets with exclusive Teclab features including reinforced shelves and backs, reinforced doors with smooth operating, easy lift-off hinges and easy to use four-point in cabinet leveling.

Yes, discovering how to put all of these heavy duty features and more into an aesthetically pleasing package was a challenge and a test of Teclab's commitment to building the best. But even we're surprised at how handsome the finished product is. In fact, a lot of our customers are saying what we've been thinking: "Why hide a cabinet that looks this good in the back-room?" Teclab has created a storage cabinet to be the leader in performance and value.

**SC-4824-GD****SIZES**

30" wide x 78" high x 18" deep (part# SC-3018*)

36" wide x 78" high x 18" deep (part# SC-3618*)

48" wide x 78" high x 18" deep (part# SC-4818*)




30" wide x 78" high x 24" deep (part# SC-3024*)

36" wide x 78" high x 24" deep (part# SC-3624*)

48" wide x 78" high x 24" deep (part# SC-4824*)

*Options (add to part number)

- WR (Clothes Bar)
- XS (Additional Shelf, 4 shelves come standard)
- GD (Glass Doors, as shown)

**SC-3018****SC-3618****SC-4818**

Teclab Wall Cabinets have always been the professional's wall cabinet and when you compare us to the competition you'll see why Teclab has been the industry leader for more than 30 years.

Designed with our exclusive Double Wall construction using heavy-duty 18 gauge, furniture grade, cold rolled steel on all parts including shelves and doors.

Tops and bottoms of all cabinets feature solid construction for added strength and durability and use of our exclusive Teclab heavy-duty wall mounting bracket system allows for simple yet solid installation.

Sturdy reinforced shelves are adjustable in 2" increments. Wall Cabinets come standard with 2 shelves, however additional shelves are available to match your storage needs.


All Teclab Wall Cabinets are available with solid double wall steel or tempered safety glass doors. Glass doors can be ordered in clear or smoked glass. We attach these doors with durable industrial grade hinges for years of trouble free use and utilize stylish stainless steel handles for a professional look.

Environmentally responsible textured powder coatings are used for a durable, long-lasting finish. Standard and special order colors allow you to personalize your laboratory.


Teclab Wall Cabinets are also offered in stainless steel for the ultimate in durability and performance.



All Teclab Wall Cabinets are available with solid steel or glass doors. Sliding glass doors are also available



TYPICAL WALL CABINET DIMENSIONS



48" wide wall cabinet

WC-4830-HSD (hinged solid doors)
WC-4830-HGD (hinged glass doors)
WC-4830-SGD (sliding glass doors)
WC-4830-SGD-PT (pass-through sliding glass doors)

42" wide wall cabinet

WC-4230-HSD (hinged solid doors)
WC-4230-HGD (hinged glass doors)
WC-4230-SGD (sliding glass doors)
WC-4230-SGD-PT (pass-through sliding glass doors)

36" wide wall cabinet

WC-3630-HSD (hinged solid doors)
WC-3630-HGD (hinged glass doors)
WC-3630-SGD (sliding glass doors)
WC-3630-SGD-PT (pass-through sliding glass doors)

30" wide wall cabinet

WC-3030-HSD (hinged solid doors)
WC-3030-HGD (hinged glass doors)
WC-3030-SGD (sliding glass doors)
WC-3030-SGD-PT (pass-through sliding glass doors)

24" wide wall cabinet

WC-2430-HSD (hinged solid doors)
WC-2430-HGD (hinged glass doors)
WC-2430-SGD (sliding glass doors)
WC-2430-SGD-PT (pass-through sliding glass doors)

18" wide wall cabinet

WC-1830-HSD-R (solid door RH pull)
WC-1830-HSD-L (solid door LH pull)
WC-1830-HGD-R (glass door RH pull)
WC-1830-HGD-L (glass door LH pull)

-XS (additional shelf)

Teclab Flipper Door Cabinets. The professional's cabinets!

That's exactly what you'll be thinking when you take a good look at the new Teclab flipper door cabinet. The Teclab design team has redefined the storage cabinet with Teclab's industry leading standard of quality and durability.

The cabinet features beefier construction with heavier material than competitive cabinets with exclusive Teclab features including reinforced sides and backs, reinforced doors with smooth operation.


Teclab has designed a flipper door cabinet to be the leader in performance and value.

**STANDARD FEATURES**

- Manufactured of heavy-duty steel (18-gauge, furniture-grade, cold-rolled)
- Solid top and bottom
- Durable textured powdercoat finish
- Extraordinary Teclab quality
- Teclab 10-year warranty
- Made in USA
- Locks and bottom shelf standard

SIZES


- 30"** wide x 16" high x 13" deep (part# FDC-30)
- 36"** wide x 16" high x 13" deep (part# FDC-36)
- 42"** wide x 16" high x 13" deep (part# FDC-42)
- 48"** wide x 16" high x 13" deep (part# FDC-48)



The Teclab Interlocking Storage cabinet. Setting a new standard in storage!

That's exactly what you'll be thinking when you take a good look at the **Teclab Interlocking Storage Cabinet**. The Teclab design team has redefined the storage cabinet with Teclab's industry leading standard of quality and durability.

The **Interlocking Storage Cabinet** features beefier construction with heavier material than competitive cabinets. With exclusive Teclab features including rugged internal frames, reinforced drawers with smooth operating slide mechanisms, easy to use four-point in-cabinet leveling and our **Teclab 10-year warranty**.

**Powdercoat Finish (Putty)****Stainless Steel****Laboratory Style Handles****Smooth strong drawers****Rear-view detail of drawers
(back removed)****Optional mobile casters accessory**

4800 SERIES, 30" or 36" wide x 48.75" high x 28" deep



C-4801-30 or 36



C-4802-30 or 36



C-4803-30 or 36



C-4804-30 or 36



C-4805-30 or 36



C-4806-30 or 36



C-4807-30 or 36



C-4808-30 or 36

6000 SERIES, 30" or 36" wide x 61.15" high x 28" deep



C-6001-30 or 36



C-6002-30 or 36



C-6003-30 or 36




C-6004-30 or 36



C-6005-30 or 36



C-6005-30 or 36





Highest quality. Best built. Most durable. Most versatile. Take your pick. They're all descriptions that have become synonymous with Teclab in the last 30 years. But our fanatical commitment to building the best is more than just talk.

The extraordinary fit and finish of Teclab Cabinet Mounted Shelving reflects that devotion to providing what customers want and need. Customers see our quality in clean, functional lines. They feel it in the solidity of a well made unit. And they experience it in the attention to the smallest of details of design assistance, installation and after-sale service.

Standard horizontal steel shelves have a stiffener formed downward to create a fully flat surface. The stiffener can also be formed upward to create a front shelf retainer.

DIMENSIONS:



		30" wide	36" wide	42" wide	48" wide
36" high complete unit	single-sided (2 shelves)	SUS-3036-2	SUS-3636-2	SUS-4236-2	SUS-4836-2
	double-sided (4 shelves)	SUS-3036-4	SUS-3636-4	SUS-4236-4	SUS-4836-4
48" high complete unit	single-sided (2 shelves)	SUS-3048-2	SUS-3648-2	SUS-4248-2	SUS-4848-2
	double-sided (4 shelves)	SUS-3048-4	SUS-3648-4	SUS-4248-4	SUS-4848-4
36" high complete unit with Electric Kit	single-sided (2 shelves)	SUS-3036-2-EP	SUS-3636-2-EP	SUS-4236-2-EP	SUS-4836-2-EP
	double-sided (4 shelves)	SUS-3036-4-EP	SUS-3636-4-EP	SUS-4236-4-EP	SUS-4836-4-EP
48" high complete unit with Electric Kit	single-sided (2 shelves)	SUS-3048-2-EP	SUS-3648-2-EP	SUS-4248-2-EP	SUS-4848-2-EP
	double-sided (4 shelves)	SUS-3048-4-EP	SUS-3648-4-EP	SUS-4248-4-EP	SUS-4848-4-EP

WALL MOUNTED SHELVING (WUS)

Standard shelf depth is 12", but other depths are available as well.

Length of shelves is variable between 18" and 96". For lengths over 48", a third(center) upright is added.

Height of Shelving Uprights is variable between 18" and 96" in 6" increments.

Shelf material can be selected from any of our Work Surface options.




	30" wide	36" wide	42" wide	48" wide
36" high complete unit	WUS-3630	WUS-3636	WUS-3642	WUS-3648
48" high complete unit	WUS-4830	WUS-4836	WUS-4842	WUS-4848

WALL MOUNTED SHELVING (WS)

Wall Standards - Used for wall mounting our steel shelves.



Shown are two 36" high Wall Standards(WS-36-W) and three 48" wide shelves.



Height	18"	24"	30"	36"	48"	72"
Part #	WS-18-W	WS-24-W	WS-30-W	WS-36-W	WS-48-W	WS-72-W

Teclab's durable floor standing shelving units are available in three widths and two depths. All are 72" high, and have 5 shelves.

Shelf height is adjustable in 3" increments.



	30" wide	36" wide	48" wide
18" deep	SU-3018-72	SU-3618-72	SU-4818-72
24" deep	SU-3024-72	SU-3624-72	SU-4824-72

Teclab has a select line of laboratory seating built to the highest standards in the industry. All seating is ergonomically designed for outstanding comfort, total body support and correct seating posture. All things leading to increased worker productivity. Each chair can be configured to your unique laboratory environment allowing you to specify such items as seating materials, chair heights and arm rests. Virtually anything to make the most comfortable seating possible.

Contact us and we'll discuss all the options available. We can help you make your next chair selection a vital addition to your laboratory while also making the most economic sense.

**PREMIUM HIGH RING ESD CHAIR & CLASS 100 CLEAN-ROOM CHAIRS**

These chairs come with a synchronized Fluid-Motion mechanism, which allows the seat and back to move in unison, providing continuous lumbar support to the user. Backrest movement "free floats" or can be locked at any point in the range. A sensitive tension control allows easy adjustment of the backrest pressure. Integral lumbar support adjusts to ensure correct seated posture. Comes standard with a 20" height adjustable foot ring, composite 5-star base and rugged 2" dual-wheel ESD casters.

**DESK HIGH CHAIR**

This chair comes standard with a seat and back that adjust independently. Backrest angle follows the user's movement by pneumatic cylinder, and can be locked at any point in the range. Integral lumbar support adjusts to ensure correct posture. Comes standard with a composite 5-star base and rugged 2" dual-wheel casters.

**PREMIUM QUALITY HIGH RING CHAIR**

This chair comes standard with a seat and back that adjust independently. Backrest angle follows the user's movement by pneumatic cylinder, and can be locked at any point in the range. Integral lumbar support adjusts to ensure correct posture. Comes standard with a 20" height adjustable foot ring, composite 5-star base and rugged 2" dual-wheel casters.

**SEATWISE HIGH RING CHAIR**

Seatwise intermediate and tall height chairs are equipped with a wide backrest that supports the user's back as they move. This feature is equipped with a tension adjustment and variable back lock. Fingertip adjustment of the seat height and backrest height adjustment with integral lumbar support make this an extremely comfortable, affordable ergonomic chair. Comes standard with a 20" diameter height adjustable foot ring, composite 5-star base and rugged 2" dual-wheel casters.

**MODERN MESH DESK CHAIR**

A contemporary mesh chair designed to create the ultimate experience in comfort, support, and durability. The modern design provides features and functionality for today's work environment.

Mesh Chair, Contoured Backrest With Lumbar Support, Adjustable Arms, Black Nylon Base, All-Purpose Casters, 18.75" - 22.25" Seat Height Adj., 275lb Weight Capacity, 5yr Warranty.

**LABORATORY STOOL**

Choose our Laboratory Stool for top performance in a wide variety of environments including: industry, education, laboratory, and health care. This stool comes standard with a 5-legged tubular base/foot ring and casters for hard floors.

**PREMIUM QUALITY HIGH RING POLYURETHANE CHAIR**

This high ring chair has durable, contoured polyurethane seat and back with ventilated ribs for lasting comfort and support. Pneumatic seat height adjustment, backrest angle and height adjustment are standard. Chair comes equipped with dual wheel rubber casters for hard floors or plastic casters for carpeted areas. Sturdy five-star reinforced plastic base with an 18" diameter adjustable chrome foot ring also standard.

TECLAB WORK SURFACES FIT VIRTUALLY ANY REQUIREMENT.

Teclab provides you with a versatile selection of work surface materials to meet any application. These materials, used in countertops, shelving and back guards, are described below.



Hardwood Butcher Block solidly constructed of edge-grained hardwood strips with clear surfaces and ends. Seasoned, kiln-dried strips are bonded together with tough resin glue and electronically cured, planed and sanded to specifications, then varnished – baked between coats – for a highly water and abrasion resistant clear finish and a satin-smooth resilient work surface. Hardwood countertops are 1 ¾" thick, shelves and back guards are 1 ¼" thick.



Teclab **Self-edged Laminate** work surfaces feature a full 1/16" high-pressure plastic laminate covering a 1 1/8" industrial grade particle board core. The bottom is coated with a phenolic backer sheet. Work surfaces and shelves are available in both 1 ¼" and 1 ¾" thickness; back guards are 1 ¼" thick. Teclab self-edged work surfaces are manufactured with matching laminate on each of the 4 edges.



Teclab **PVC-edged Laminate** work surfaces are constructed in the same manner as Self-Edged Laminate work surfaces but with a 3mm thick PVC edge banding to provide excellent resistance to bumps and scrapes that could damage the edge laminate and cause it to peel away. The edge-banding color closely matches the laminate color and comes standard on Teclab Work Stations.



Teclab **Composite Resinwood** work surfaces are constructed with a 1" thick industrial grade particle board core. The top and bottom are 3/8" thick urea formaldehyde-free phenolic resinwood with a two-coat catalyzed varnish finish. The front edge has a 3/8" radius. Overall thickness is 1 ¾".



Chem-Surf® Laminate work surfaces are used in environments where the possibility of chemical spills rules out the use of standard laminates on work surfaces or where costs prohibit epoxy resin, slate, marble or granite surfaces.



Teclab **Electrostatic Dissipative (ESD)** work surfaces are manufactured using ESD laminates. A 180° roll on the front edge ensures continuous conductivity along the front edge of the work surface. These countertops provide a colorful, durable, and reliable work surface with the recommended surface resistivity of 1×10^6 to $1 \times 10^9 \Omega$.



Tec-Resin Epoxy work surfaces are monolithic, offering consistent chemical resistance throughout. They are molded from modified epoxy resin, compounded and heat cured. Tec-Resin Epoxy offers outstanding resistance to acid, alkali and solvents, as well as physical abuse. The low-sheen black satin finish makes Tec-Resin Epoxy an ideal work surface, both aesthetically and functionally. Available in Black, Grey, or White.



Teclab **Phenolic Resin** Work Surfaces provide an ideal choice in environments where hygiene, avoidance of contamination, sustainability, ergonomics and safety are of importance

Available in three colors: White, Gray and Black.



Teclab **Stainless Steel** work surfaces are constructed of 304 stainless steel with a #4 finish. Corners are formed, welded, then polished and ground to match the pattern. Reinforcing channel members are welded to the underside to provide structural rigidity to the entire surface. Overall thickness is 1".


WL-1800

Bench-Height Base Cabinets

18" wide x 35.5" high x 22" deep



Nearly any base cabinet can be made **MOBILE**.



WL-1801



WL-1802



WL-1803



WL-1804



WL-1805



WL-1806



WL-1807



WL-1808



WL-1809



WL-1810



WL-1811



WL-1812

WL-1800 Bench-Height Base Cabinets, Cont.



WL-1813



WL-1814



WL-1815



WL-1816



WL-1817



WL-1818



WL-1819



WL-1820



WL-1821




WL-1822



WL-1823



WL-1824



WL-1826



WL-2400

Bench-Height Base Cabinets


24" wide x 35.5" high x 22" deep




Nearly any base cabinet can be made **MOBILE**.



WL-2401




WL-2402




WL-2403



WL-2404



WL-2405



WL-2406




WL-2407



WL-2408



WL-2409



WL-2410



WL-2411



WL-2412

WL-2400 Bench-Height Base Cabinets, Cont.



WL-2413



WL-2414



WL-2415



WL-2416



WL-2417



WL-2418



WL-2419



WL-2420



WL-2421




WL-2422



WL-2423



WL-2424



WL-2426



WL-3000

Bench-Height Base Cabinets

30" wide x 35.5" high x 22" deep



Nearly any base cabinet can be made **MOBILE**.



WL-3001




WL-3002



WL-3003



WL-3004



WL-3005




WL-3006




WL-3007




WL-3008



WL-3009




WL-3010




WL-3011



WL-3012

WL-3000 Bench-Height Base Cabinets, Cont.

WL-3014



WL-3015



WL-3016



WL-3017



WL-3018


WL-3600

Bench-Height Base Cabinets

36" wide x 35.5" high x 22" deep



Nearly any base cabinet can be made **MOBILE**.



WL-3601



WL-3602



WL-3603



WL-3604



WL-3605



WL-3606



WL-3607



WL-3608



WL-3609



WL-3610



WL-3611



WL-3612

WL-3600 Bench-Height Base Cabinets, Cont.



WL-3613



WL-3614



WL-3615



WL-3616




WL-3617



WL-3618



WL-3619



WL-3620



WL-3621



WL-3622



WL-3623



WL-3624



WL-3625



WL-3626



WL-3627



WL-3628

WL-3600 Bench-Height Base Cabinets, Cont.



WL-3629



WL-3630



WL-3631



WL-3632



WL-3633



WL-3634



WL-3635



WL-3636



WL-3637



WL-3638



WL-3639



WL-3640



WL-3641



WL-3642



WL-3643



WL-3644

WL-3600 Bench-Height Base Cabinets, Cont.



WL-3645



WL-3646



WL-3647



WL-3648



WL-3649



WL-3650



WL-3651



WL-3652



WL-3653



WL-3654



WL-3655



WL-3656


WL-4200

Bench-Height Base Cabinets

42" wide x 35.5" high x 22" deep



Nearly any base cabinet can be made **MOBILE**.



WL-4201



WL-4202



WL-4203



WL-4204



WL-4205



WL-4206



WL-4207



WL-4208



WL-4209



WL-4210



WL-4211



WL-4212

WL-4200 Bench-Height Base Cabinets, Cont.



WL-4213



WL-4214



WL-4215



WL-4216



WL-4217



WL-4218



WL-4219



WL-4220



WL-4221



WL-4222



WL-4223



WL-4224



WL-4225



WL-4226



WL-4227



WL-4228

WL-4200 Bench-Height Base Cabinets, Cont.



WL-4229



WL-4230



WL-4231



WL-4232



WL-4233



WL-4234



WL-4235



WL-4236



WL-4237



WL-4238



WL-4239



WL-4240



WL-4241



WL-4242



WL-4243



WL-4244

WL-4200 Bench-Height Base Cabinets, Cont.



WL-4245




WL-4246



WL-4247



WL-4248



WL-4249



WL-4250



WL-4251



WL-4252



WL-4253


WL-4800

Bench-Height Base Cabinets

48" wide x 35.5" high x 22" deep



Nearly any base cabinet can be made **MOBILE**.



WL-4801



WL-4802



WL-4803



WL-4804



WL-4805



WL-4806



WL-4807



WL-4808



WL-4809



WL-4810



WL-4811



WL-4812



WL-4813



WL-4814



WL-4815



WL-4816

WL-4800 Bench-Height Base Cabinets, Cont.



WL-4817



WL-4818



WL-4819



WL-4820



WL-4821



WL-4822



WL-4823



WL-4824



WL-4825



WL-4826



WL-4827



WL-4828



WL-4829



WL-4830



WL-4831



WL-4832



WL-4833



WL-4834



WL-4835



WL-4836

WL-4800 Bench-Height Base Cabinets, Cont.



WL-4837



WL-4838



WL-4839



WL-4840



WL-4841



WL-4842



WL-4843



WL-4844



WL-4845




WL-4846



WL-4847



WL-4848



WL-4849



WL-4850



WL-4851



WL-4852L



WL-4852R




WL-4853

END ISLAND

Bench-Height Base Cabinets

Various widths x 35.5" high x 22" deep



WL-4602



WL-4603



WL-4604



WL-4605



WL-5802



WL-5819




WL-4601
SINK CABINET




WL-5201
SINK CABINET




WL-5801
SINK CABINET



WL-6101
SINK CABINET



WL-6401
SINK CABINET



WL-7001
SINK CABINET

DL-1800


Desk-Height Base Cabinets

18" wide x 28" high x 22" deep



Nearly any base cabinet can be made **MOBILE**.*

*Increases the height of Desk-High cabinets by 3".



DL-1801



DL-1802



DL-1803



DL-1804



DL-1805



DL-1806



DL-1807



DL-1808



DL-1809



DL-1810



DL-1811



DL-1812



DL-1813



DL-1814



DL-1816



DL-1817



DL-1818

DL-2400



Desk-Height Base Cabinets

24" wide x 28" high x 22" deep




Nearly any base cabinet can be made **MOBILE**.*

*Increases the height of Desk-High cabinets by 3".



DL-2401



DL-2402



DL-2403



DL-2404



DL-2405



DL-2406



DL-2407



DL-2408



DL-2409



DL-2410



DL-2411



DL-2412



DL-2413



DL-2414



DL-2416



DL-2417



DL-2418



DL-2419

DL-3000


Desk-Height Base Cabinets

30" wide x 28" high x 22" deep



Nearly any base cabinet can be made **MOBILE**.*

*Increases the height of Desk-High cabinets by 3".



DL-3001



DL-3002



DL-3003



DL-3004



DL-3005



DL-3006



DL-3007



DL-3008



DL-3010

DL-3600


Desk-Height Base Cabinets

36" wide x 28" high x 22" deep



Nearly any base cabinet can be made **MOBILE**.*

*Increases the height of Desk-High cabinets by 3".



DL-3601



DL-3602



DL-3603



DL-3604



DL-3605



DL-3607




DL-3608




DL-3609

HANGING CABINETS

18" wide x various height x 22" deep, Hanging Cabinets



WUTD-1800




WUTS-1801



WUTS-1802




WUTS-1803



WUTS-1804



WUTS-1805



WUTS-1806-RHP



WUTS-1806-LHP



WUTS-1807



WUTS-1808



WUTS-1809

ACID STORAGE CABINETS

Acid Storage Cabinets safely store and vent acids and other corrosive liquids.

ASC-18-RH (18" wide Acid Storage Cabinet with right-hand hinge.)

ASC-18-LH (18" wide Acid Storage Cabinet with left-hand hinge.)


ASC-24-RH (24" wide Acid Storage Cabinet with right-hand hinge.)

ASC-24-LH (24" wide Acid Storage Cabinet with left-hand hinge.)

ASC-30 (30" wide Acid Storage Cabinet with double doors.)

ASC-36 (36" wide Acid Storage Cabinet with double doors.)

ASC-48 (48" wide Acid Storage Cabinet with double doors.)

**SOLVENT STORAGE CABINETS**

Solvent Storage Cabinets safely store and vent solvents and other flammable liquids. In case of a building fire, these cabinets protect the contents inside them.


FSC-24-RH (24" wide Solvent Storage Cabinet with right-hand hinge.)

FSC-24-LH (24" wide Solvent Storage Cabinet with left-hand hinge.)

FSC-30 (30" wide Solvent Storage Cabinet with double doors.)

FSC-36 (36" wide Solvent Storage Cabinet with double doors.)

FSC-48 (48" wide Solvent Storage Cabinet with double doors.)




CORNER ADAPTERS

Corner Adapters provide a flush filler panel when placed between Corner Cabinet and adjacent Base Cabinet or Panel Leg. Mounts directly to cabinet side for easy installation and is finished in the same manner as the Base Cabinets.

WCA-1500 (Bench Height (35 ½" high x variable width)


WCA-1000 (Desk Height (28" high x variable width)

**FILLER PANELS**

Filler Panels provide a smooth filler between the cabinet and wall. Each Filler Panel is finished in the same manner as the Base Cabinets.

WFP-1500 (Bench Height (35 ½" high x variable width)


WFP-1000 (28" high x variable width)

**CHANNEL CABINET SPACERS**

Placing the Channel Cabinet Spacers between the Base Cabinets provides a flush filler panel. Each Channel Cabinet Spacer is finished in the same manner as the Base Cabinets.


WCCS-15 (Bench Height (35 ½" high x variable width)

WCCS-10 (Desk Height (28" high x variable width)

**EXTRA SHELF FOR DOOR UNITS**

Available for any Base Cabinet with a door to provide added storage space. Each Shelf is finished in the same manner as the Base Cabinet.

xx-xxxx-XS (To order, simply add -XS after Base Cabinet model number.)



MODESTY PANELS

Each Modesty Panel is a reinforced assembly equipped with a wiring channel and is finished in the same manner as the Base Cabinets. Modesty Panels extend between Base Cabinets or Panel Legs to provide stability to some workbenches. Modesty Panels also provide some privacy and can be used to hide plumbing, wiring or duct work behind the workbench.

Modesty Panels - Full Height (shown)


WMP-15xx (Bench Height (35 ½" h x 1" d x variable width)

WMP-10xx (Desk Height (28" h x 1" d x variable width)

Modesty Panels - For use with Apron Drawers

WMPA-15xx (Bench Height (29 ¾" h x 1" d x variable width)

WMPA-10xx (Desk Height (23 ¼" h x 1" d x variable width).

**APRON DRAWERS**

Apron Drawers utilize the standard "A" size Base Cabinet drawer in the six available widths. All Apron Drawers are manufactured to accept locks. Manufactured and finished to the same rigid specifications as the Base Cabinets.

WAPD-18 (Case: 18" wide x 5" high x 22" deep)


WAPD-24 (Case: 24" wide x 5" high x 22" deep)

WAPD-30 (Case: 30" wide x 5" high x 22" deep)


WAPD-36 (Case: 36" wide x 5" high x 22" deep)

WAPD-42 (Case: 42" wide x 5" high x 22" deep)

WAPD-48 (Case: 48" wide x 5" high x 22" deep)

**APRON FILLERS**

Apron Fillers are used to fill the gap between fixed-sized Apron Drawers and Base Cabinets. Manufactured and finished to the same rigid specifications as the Base Cabinets.



PANEL AND SUPPORT LEGS

Panel Legs are constructed of heavy-gauge, welded steel assembly, with easily removable side panels. Each is provided with knockouts in the back frame and punch outs in the top to allow for flexibility during installation. Each Panel Leg is finished in the same manner as the Base Cabinet.

Support Legs (not shown) are used to hold up work surfaces which require added support and are located away from walls.

Panel Leg WPL-15

(Bench Height (35 ½" high x 4" wide x 22" deep)

Panel Leg WPL-10


(28" high x 4" wide x 22" deep)

Support Leg WLEG-15

(Bench Height (35 ½" high x 4" wide x 8" deep)

Support Leg WLEG-10

(28" high x 4" wide x 8" deep)

**PANEL LEGS WITH ELECTRICAL POWER KITS**


Panel Legs are available with a variety of power supplies. Standard kit, duplex, or it can be customized to fit any power supply needed.

WPL-15-EP

(Bench Height (35 ½" high x 4" wide x 22" deep)


WPL-10-EP

(28" high x 4" wide x 22" deep)



POWERED BACKGUARDS

Powered Back Guards (variable length x 4" H x 2" D) are available in lengths to match the length of the work surface or shelf and come standard with Electrical Power Kit EP-1 (other EP kits are available). **PBG-72-EP** is 72" long providing 6 electrical outlets in EP-1 kit, **PBG-60-EP** is 60" long providing 5 electrical outlets in EP-1 kit, and **PBG-96-EP** is 96" long providing 8 electrical outlets in 2 EP-1 kits.)




VERTICAL POWER POLES

Vertical Power Poles (variable height x 7" W x 3" D) are used to contain utilities originating from the ceiling such as electric, plumbing or gas.

VPP-60 (60" Height)

VPP-96 (96" Height)




COUNTERTOP SUPPORT BRACKETS

Countertop Support Brackets are used to support work surfaces along a wall or to connect multiple work surfaces. They are built from heavy-gauge steel and supplied with the necessary attaching hardware.

CSB-18 (18" long x 4" deep)

CSB-18-90 (18" long, 2" x 2" L shaped)

CSB-18-90-W (18" long, 2" x 2" L shaped + adjustable levelers)



CASTERS

Most Teclab Base Cabinets can be used as Mobile Base Cabinets with the addition of Casters. Standard Casters are 4" high (3" wheels).

The height of Laboratory style Mobile bench-height cabinets match the height of fixed cabinets. The height of Laboratory style Mobile desk-height cabinets increases by 3" inches.




C-1500

Bench-Height Base Cabinets

15" wide x 33.75" high x 28" deep



Nearly any base cabinet can be made **MOBILE**.



C-1501



C-1502



C-1503




C-1504




C-1505



C-1506



C-1507



C-1508



C-1509



C-1510



C-1511



C-1512



C-1513



C-1514



C-1515



C-1516



C-1517



C-1518

C-1500 Bench-Height Base Cabinets, Cont.



C-1519



C-1520



C-1521



C-1522



C-1523



C-1524



C-1525



C-1526



C-1527



C-1528



C-1529



C-1530



C-1531



C-1532



C-1533



C-1534



C-1535



C-1536



C-1537



C-1538



C-1539



C-1540



C-1541



C-1542


C-2200

Bench-Height Base Cabinets

22" wide x 33.75" high x 28" deep



Nearly any base cabinet can be made **MOBILE**.



C-2201



C-2202



C-2203



C-2204



C-2205



C-2206



C-2207



C-2208



C-2209



C-2210



C-2211



C-2212

C-2200 Bench-Height Base Cabinets, Cont.



C-2213



C-2214



C-2215



C-2216



C-2217



C-2218



C-2219



C-2220



C-2221



C-2222



C-2223



C-2224



C-2225



C-2226



C-2227



C-2228


C-3000

Bench-Height Base Cabinets


30" wide x 33.75" high x 28" deep



Nearly any base cabinet can be made **MOBILE**.



C-3001



C-3002



C-3003



C-3004



C-3005



C-3006



C-3007



C-3008



C-3009



C-3010



C-3011



C-3012

C-3000 Bench-Height Base Cabinets, Cont.



C-3013



C-3014



C-3015



C-3016



C-3017



C-3018



C-3019



C-3020



C-3025


C-3600

Bench-Height Base Cabinets

36" wide x 33.75" high x 28" deep



Nearly any base cabinet can be made **MOBILE**.



C-3601



C-3602



C-3603



C-3604



C-3605



C-3606



C-3607



C-3608



C-3609



C-3610



C-3611



C-3612



C-3613



C-3614

C-1000


Bench-Height Base Cabinets

15" wide x 27.75" high x 28" deep



Nearly any base cabinet can be made **MOBILE**.*

*Increases the height of Desk-High cabinets by 2 3/4".



C-1001



C-1002



C-1003



C-1004



C-1005



C-1006



C-1007



C-1008



C-1009



C-1010



C-1011



C-1012




C-1013



C-1014



C-1015



C-1016



C-1017



C-1018



C-1019

C-2100


Desk-Height Base Cabinets

22" wide x 27.75" high x 28" deep



Nearly any base cabinet can be made **MOBILE**.*

*Increases the height of Desk-High cabinets by 2 3/4".



C-2101



C-2102



C-2103



C-2104



C-2105



C-2106



C-2107



C-2108



C-2109



C-2110



C-2111



C-2112



C-2113



C-2114



C-2115



C-2116



C-2117

C-2000


Bench-Height Base Cabinets

30" wide x 33.75" high x 28" deep



Nearly any base cabinet can be made **MOBILE**.*

*Increases the height of Desk-High cabinets by 2 3/4".



C-2001



C-2002



C-2003



C-2004



C-2005



C-2006



C-2007



C-2008



C-2009



C-2010



C-2011



C-2012



C-2013



C-2014

C-3700


Bench-Height Base Cabinets

36" wide x 33.75" high x 28" deep



Nearly any base cabinet can be made **MOBILE**.*

*Increases the height of Desk-High cabinets by 2 3/4".



C-3701



C-3702



C-3703



C-3704



C-3705



C-3706



C-3707



C-3708



C-3709


C-6000

Stacked Base Cabinets

15", 22" or 30" wide x 61.5" high x 28" deep



Nearly any base cabinet can be made **MOBILE**.



C-6001-15



C-6002-15



C-6003-15



C-6004-15



C-6005-15



C-6006-15



C-6007-22



C-6008-22



C-6009-22



C-6010-22



C-6011-22



C-6012-22



C-6013-30



C-6014-30



C-6015-30



C-6016-30



C-6017-30




C-6018-30

DC-1500

Double-sided Island Base Cabinets


15" wide x 33.75" high x 34, 40 or 46" deep



Nearly any base cabinet can be made **MOBILE**.



DC-1501-36
DC-1501-42
DC-1501-48



DC-1502-36
DC-1502-42
DC-1502-48



DC-1506-36
DC-1506-42
DC-1506-48



DC-1509-36
DC-1509-42
DC-1509-48



DC-1522-36
DC-1522-42
DC-1522-48



DC-1531-36
DC-1531-42
DC-1531-48



DC-1532-36
DC-1532-42
DC-1532-48



DC-1534-36
DC-1534-42
DC-1534-48



DC-1535-36
DC-1535-42
DC-1535-48


DC-2200

Double-sided Island Base Cabinets

22" wide x 33.75" high x 34, 40 or 46" deep




Nearly any base cabinet can be made **MOBILE**.



DC-2203-36
DC-2203-42
DC-2203-48



DC-2204-36
DC-2204-42
DC-2204-48



DC-2210-36
DC-2210-42
DC-2210-48



DC-2214-36
DC-2214-42
DC-2214-48




DC-2218-36
DC-2218-42
DC-2218-48



DC-2219-36
DC-2219-42
DC-2219-48



DC-2221-36
DC-2221-42
DC-2221-48



DC-2222-36
DC-2222-42
DC-2222-48


DC-3000

Double-sided Island Base Cabinets

30" wide x 33.75" high x 34, 40 or 46" deep



Nearly any base cabinet can be made **MOBILE**.



DC-3001-36
DC-3001-42
DC-3001-48



DC-3012-36
DC-3012-42
DC-3012-48



DC-3013-36
DC-3013-42
DC-3013-48



DC-3015-36
DC-3015-42
DC-3015-48



DC-3016-36
DC-3016-42
DC-3016-48





DC-3017-36
DC-3017-42
DC-3017-48



DC-3019-36
DC-3019-42
DC-3019-48

HANGING CABINETS

15" wide x various height x 28" deep



UTD-1500




UTS-1501




UTS-1502



UTS-1503



UTS-1504



UTS-1505



UTS-1506-RHP



UTS-1506-LHP



UTS-1507



UTS-1508




UTS-1510

CORNER ADAPTERS

Corner Adapters provide a flush filler panel when placed between Corner Cabinet and adjacent Base Cabinet or Panel Leg. Mounts directly to cabinet side for easy installation and is finished in the same manner as the Base Cabinets.

CA-1500 (Bench Height (33 ¾" high x variable width)


CA-1000 (Desk Height (27 ¾" high x variable width)

**FILLER PANELS**

Filler Panels provide a smooth filler between the cabinet and wall. Each Filler Panel is finished in the same manner as the Base Cabinets.

FP-1500 (Bench Height (33 ¾" high x variable width)


FP-1000 (27 ¾" high x variable width)

**CHANNEL CABINET SPACERS**

Placing the Channel Cabinet Spacers between the Base Cabinets provides a flush filler panel. Each Channel Cabinet Spacer is finished in the same manner as the Base Cabinets.


CCS-15 (Bench Height (33 ¾" high x variable width)

CCS-10 (Desk Height (27 ¾" high x variable width)

**EXTRA SHELF FOR DOOR UNITS**

Available for any Base Cabinet with a door to provide added storage space. Each Shelf is finished in the same manner as the Base Cabinet.

xx-xxxx-XS (To order, simply add -XS after Base Cabinet model number.)



APRON DRAWERS


Apron Drawers are available in four widths:

APD-1500 is 15" W x 3 7/32" H x 27" D


APD-2200 is 22" W x 3 7/32" H x 27" D

APD-3000 is 30" W x 3 7/32" H x 27" D


APD-3600 is 36" W x 3 7/32" H x 27" D)

**APRON FILLERS**

Apron Fillers are used to fill the gap between fixed-sized Apron Drawers and Base Cabinets. Manufactured and finished to the same rigid specifications as the Base Cabinets.

**SHELF SUPPORTS**

Shelf Supports are available in many heights, with or without Electrical Power Kit EP-1 (other EP kits are available). These rigid supports are constructed of a heavy-gauge welded steel assembly, with easily removable side panels that provide easy access to wiring. Each Shelf Support (EP) has knockouts in the back and punch-outs in the top and bottom to allow for flexibility during installation. Finished in the same manner as the Base Cabinets.



SS-12 (12" H x 4" W x 8" D)

SS-12-EP (12" H x 4" W x 8" D)

SS-16 (16" H x 4" W x 8" D)

SS-16-EP (16" H x 4" W x 8" D)

SS-18 (18" H x 4" W x 8" D)

SS-18-EP (18" H x 4" W x 8" D)

SS-21 (21" H x 4" W x 8" D)

SS-21-EP (21" H x 4" W x 8" D)


SS-24 (24" H x 4" W x 8" D)

SS-24-EP (24" H x 4" W x 8" D)


FOOTRESTS

Footrests (**FR6-xx**) are reinforced heavy-gauge steel and are supplied with the necessary attaching hardware. Base Cabinets and Panel Legs are pre-punched to simplify assembly. Foot Rests are finished in black.

Typical dimensions: 1" H x 6" D x variable Width.

**PANEL LEGS**

Panel Legs are available in two heights to correspond with the two available heights of standard Base Cabinets. These Panel Legs are constructed of heavy-gauge, welded steel assembly, with easily removable side panels for access to wiring. Each is provided with knockouts in the back frame and punch outs in the top to allow for flexibility during installation. Each Panel Leg (EP) is furnished with Electrical Power Kit EP-1 (other EP kits are available), and is finished in the same manner as the Base Cabinets. Should be used in conjunction with Foot Rests and Modesty Panels to add rigidity to any Work Station assembly.



PL-15 (33 1/4" H x 4" W x 28" D)

PL-10 (27 3/4" H x 4" W x 28" D)

PL-15-EP (33 1/4" H x 4" W x 28" D)


PL-10-EP (27 3/4" H x 4" W x 28" D)

SUPPORT LEGS

Support Legs are available in two heights to correspond with the two available heights of standard Base Cabinets. These Support Legs are constructed of heavy-gauge steel, and are finished in the same manner as the Base Cabinets.


LEG-15 (33 1/4" H x 4" W x 8" D)

LEG-10 (27 3/4" H x 4" W x 8" D)




POWERED BACK GUARDS

Powered Back Guards (variable length x 4" H x 2" D) are available in lengths to match the length of the work surface or shelf and come standard with Electrical Power Kit EP-1 (other EP kits are available). **PBG-72-EP** is 72" long providing 6 electrical outlets in EP-1 kit, **PBG-60-EP** is 60" long providing 5 electrical outlets in EP-1 kit, and **PBG-96-EP** is 96" long providing 8 electrical outlets in 2 EP-1 kits.

**MODESTY PANELS**

Modesty Panels (23 7/8" H x 1" D x variable width) Each Modesty Panel is a reinforced welded assembly equipped with a wiring channel and is finished in the same manner as the Base Cabinets. Modesty Panels extend between Base Cabinets or Panel Legs to provide stability to some workbenches. Modesty Panels also provide some privacy and can be used to hide plumbing, wiring or duct work behind the workbench.


**COUNTERTOP SUPPORT BRACKETS**

Countertop Support Brackets are used to support work surfaces along a wall or to connect multiple work surfaces. They are built from heavy-gauge steel and supplied with the necessary attaching hardware.

CSB-18 (18" long x 4" deep)


CSB-18-90 (18" long, 2" x 2" L shaped)

CSB-18-90-W (18" long, 2" x 2" L shaped + adjustable levelers).




DRAWER DIVIDER KITS

Organize your Teclab cabinets with these drop-in kits. Made to fit "A" and "B" drawers in all cabinet widths. Each kit comes with 2 side rails, a drawer mat, and enough spacers to configure drawers in many ways.


**HANGING FILE ADAPTER KITS**

Turn any "C" drawer into a file drawer with these helpful kits.

**UNDER-SHELF LINEAR LED LIGHT KITS**


Under Shelf Light Kits install under Shelves or Cabinets and provide diffused light to the work surface.

Lengths available: 12", 17", 31", 44" and 58"

**MAGNIFYING TASK LIGHT**


The 5", 3-diopter magnifier lens is illuminated by powerful/dim-mable LED light-source. It has heavy-duty, all-metal construction with a sleek design, adjustable arm and excellent light output making it suitable for a variety of industrial applications. The clamp-on base holds the lamp securely to the work surface. UL listed. 115/120v.

part# **LMG-2230**



LOCKS


Laboratories can contain some extremely valuable items, so for added security all drawers and doors in a Teclab casework cabinet can be fitted with locks. All locks can be keyed in various combinations.

**SMALL PARTS DRAWER CABINETS**

SP-8106 Small Parts Drawer Unit is equipped with 6 drawers.

SP-8118 Small Parts Drawer Unit is equipped with 18 drawers.


Small Parts Drawer Spacers: **XP-8104** Box of 12 each for use with SP-8118 and SP-8106 Small Parts Drawer Cabinets.

**LABEL HOLDERS**

Our aluminum label holders can be mounted on any door or drawer. This provides the most efficient means of keeping track of items in your laboratory. With a handsome brushed finish, our label holders are useful additions to every laboratory.


LH-35175 for 3" x 2" cards.

LH-55375 for 5" x 3" cards.



DRYING RACKS

Drying Racks are an important addition to any work station to keep your laboratory organized allowing you to operate more efficiently. Teclab has a full line of Stainless Steel and Epoxy Resin drying racks with built-in drip collector trays to keep work areas clean. Available in various board sizes and multiples of pegs to suit your laboratory needs. Our Stainless version is made of durable polished 304 Stainless Steel for years of trouble free use. The Epoxy Resin pegboards feature a no twist oval base design set at 30° for maximum support. Both Stainless Steel and Epoxy Resin drying racks are easy to install and clean thanks to removable pegs. Acrylic pegboards are available on request. Drying Racks also have a variety of mounting options too.



FAUCETS

Mixing Faucets:

FM8-412 Hot & Cold Mixing Faucet with 8" Gooseneck Spout

FM8-412-EW Mixing Faucet with 8" Gooseneck Spout & Eye-wash

FM8-412-MW Wall Mounted Mixing Faucet with 8" Goose-neck Spout

FM8-412-VB Vacuum Breaker Mixing Faucet


Single Faucets:

FS8-612 Cold Water Faucet with 8" Gooseneck Spout

FS8-612-DI Distilled/Deionized Water Faucet with 8" Goose-neck Spout

F6-412 Deck Mounted Single Spout (for use with foot pedal valves)

xx-xxxx-FP Foot Pedal Valves (Xs are the cabinet part#)




SAFETY WASHES


EW-10.5-SA Deck Mounted Swing Away Eye Wash

SS-EW Combo Safety Shower, Eye/Face Wash

EW-HH Deck Mounted Hand Held Eye Wash

**SINKS****Stainless Steel Sinks**

Teclab's stainless steel sinks are available in single and double bowl styles. A large choice of sink sizes makes outfitting your laboratory easy. Seamlessly formed of 18 gauge, 304 nickel bearing stainless these durable sinks are hand blended to a lustrous satin finish. To prevent condensation and dampen sound the underside is fully undercoated. All stainless steel sinks are furnished with 1 faucet hole.

**Drop-in Resin Sinks**

Teclab's Epoxy Resin drop in sink gives you design flexibility and eliminates the need for bulky under-sink supports. The weight bearing lip surrounding the sink is set into the recessed opening of the work surface providing a streamlined look. It also makes this sink compliant with the Americans with Disabilities Act. The rounded exterior edges assist in the users comfort and safety. Rear corner drains also allow for placement of plumbing to the back of the cabinet freeing up valuable under counter storage space.

GAS/AIR/VACUUM NOZZLES

Deck Mounted:

GAV-D-SN Single Nozzle

GAV-D-DN Double Nozzle

GAV-D-QN Quad Nozzle


Panel Mounted:

GAV-P-SN Single Nozzle

GAV-P-DN Double Nozzle


Wall Mounted:

GAV-W-SN Single Nozzle



WORK SURFACE MOUNTED ELECTRICAL PEDESTAL ASSEMBLIES

Electrical Pedestal Assemblies are mounted on work surfaces virtually anywhere you need them making all your electrical connections easier. Switching from tool to tool is made effortless without the need to bend down to floor level just to change power sources. Available in several combinations of outlets and voltages, the pedestal assemblies are also UL listed and have provisions for a grounding screw. Each electrical outlet box has a brushed aluminum finish and 3/4" NPSM mounting shank assembly. Keep these items in mind when planning your next laboratory.



ELECTRICAL POWER KITS


EP-1 Kit includes 4, 5, or 6 15 amp, 110 volt grounded outlets and illuminated combination circuit breaker rocker switch, 9ft. power cord, plus all necessary components for assembly.

EP-2 Kit includes 15 amp duplex outlets and illuminated combination circuit breaker rocker switch, 9ft. power cord, plus all necessary components for assembly.

EP-3 Kit includes 20 amp, 110 volt duplex outlets.


EP-4 Kit includes 15 amp, 220 volt duplex outlets.

EP-GFI Kit includes Ground Fault Interrupter (GFI) outlets.



KEYBOARD TRAY

Our Keyboard Tray (**KBT-30-MT**) has tilt adjustment and is mounted on a slide. It includes a movable mouse tray suitable for left or right handed use.

**MONITOR ARMS**

Monitor Arms come in several variations to suit your needs.

MAHD-3500 - 18" Radial Arm Monitor Holder

MAHD-7500 - 27" Radial Arm Monitor Holder


MAHD-7509 - 27" Radial Arm Monitor w/ Keyboard Tray

MAHD-9120 - Side by Side Monitor Holder

MAHD-9120-Q - Quad Monitor Holder

MAHD-1000 - Fixed Monitor Holder


(For TES Bench Cross-Supports)

**CPU HOLDERS**

CPU Holders are 20" high and have adjustable width.

CPU-1020 (Mounts to underside of work surface and to side of base cabinet.)

CPU-1020-M (Both halves mount to underside of work surface.)



POWDER COATED FINISH


Teclab furniture is available in environmentally responsible textured powder coatings for a durable, long-lasting finish. Standard, special order, and even custom paint matching allow you to personalize your order.



* Denotes a non-standard color. ** Custom paint color matching is available


BRUSHED STAINLESS STEEL

Teclab cabinets and many other of our products are available constructed of stainless steel. Stainless steel is strong, long-lasting and easy to sanitize. It's resistant to heat and corrosion, mold, fungi and bacteria.



WORK SURFACE COLOR OPTIONS

Visual comparison. Color choices vary depending on work surface type.



NOTE: Colors in this document are approximate. Contact us for color samples.