

Teclab Preinstallation Checklist.

Thank you for choosing Teclab factory installation with your order. There is some general information, and a few questions, that we have about your specific environment which will help our installers complete the job on-time and with minimal disruptions to your company. Please answer all questions below. Write "NA" for anything that does not apply, and feel free to contact us with any questions or clarifications.

Customer Information

Name and contact information of person completing form:

First Name*	Last Name*
Company Name*	Address
City	State
Zip Code	
Phone*	Email Address*

General Information

Reporting Location:

PPE Requirements:

COVID Restrictions:

Special Security Requirements:

Work Hours and Days Available:

Special Situations, Conditions:

If possible, please provide photos of the area the furniture is to be installed. Check if photos will be provided.

☐

Please verify the condition, and location of the furniture. In order to start the installation we will need verification that the furniture has been received in good order and has been placed in, or directly outside of the installation area.

Please verify that no modifications have been made to the installation area since the plans were finalized.

Will there be any other trades working in the same area? If yes, what trades and what will they be working on?

What is the material of the walls where the installation will take place?

Are there any obstacles on the walls we should be aware of? If yes, please explain.

Electrical Information

Teclab provides one 8' power cord per powered component. Do you have any specific needs for routing or concealing power cords?

What is the height to the center of the outlets on the wall of the installation area?

Since each powered component has its own cord, there can be multiple power cords per bench. Do you need the number of power cords reduced per bench? If so, please explain the desired configuration.

Plumbing Information

Teclab sink cabinets have a removable back to allow for entry and exit of water and drain lines.

Will the plumbing connections be entering the cabinet from the wall or the floor?

What are the rough dimensions on the wall for the plumbing stubwork?

Are there any special circumstances that may require cabinet modification for plumbing access?

Additional Information

- ☐ High Security Building
- ☐ Residential
- ☐ Elevator
- ☐ Cert Of Insurance Property Mgmt
- ☐ Cert Of Insurance End User
- ☐ Over 50 Yard Push
- ☐ Street Offload
- ☐ Union Building
- ☐ Floor Protection
- ☐ Wall Protection Elevator
- ☐ Cert Of Insurance Landlord
- ☐ Permits Required
- ☐ No Loading Dock

of Steps (12 steps = 1 flight)

of Flights