


BENCH HIGH BASE CABINET C-1505

Teclab Base Cabinets are constructed of only top quality, heavy-duty, furniture-grade, cold-rolled steel, double-braced and welded at all critical joints; twin full-perimeter frames formed of heavy gauge steel, welded to provide strength and rigidity; heavy-duty drawers that glide effortlessly on precision steel ball bearings; cabinet doors that open and close with authority and precision using self-closing hinges; brushed aluminum pulls that are recessed for clean lines and unobtrusive operation; durable, textured paint finish, in a wide selection of mar-resistant, easy-to-clean colors.

Any Base Cabinet can be made **MOBILE**.



TOP

C-1505 FEATURES

One "A" door with RH pull and sliding bottom tray, one "C" drawer base cabinet.

Cabinet Dimensions: 15" wide, 33 3/4" high, 28" deep

"A" Door with Right Hand Pull

Door Face Dimensions: 15" wide, 13 1/2" high, 1" thick.

Sliding Tray

Dimensions: 12" wide, 2" high, 25 1/2" deep.

Weight Capacity: 75 lbs.


"C" Drawer

Drawer Face Dimensions: 15" wide, 13 1/2" high, 1" thick.

Drawer Body Inside Dimensions: 12" wide, 9 1/2" high, 25 1/2" deep.

Weight Capacity: 75 lbs.

All 15" wide "C" drawers accept standard hanging files.



FRONT



SIDE