


BENCH HIGH BASE CABINET C-1520-EP-1

Teclab Base Cabinets are constructed of only top quality, heavy-duty, furniture-grade, cold-rolled steel, double-braced and welded at all critical joints; twin full-perimeter frames formed of heavy gauge steel, welded to provide strength and rigidity; heavy-duty drawers that glide effortlessly on precision steel ball bearings; cabinet doors that open and close with authority and precision using self-closing hinges; brushed aluminum pulls that are recessed for clean lines and unobtrusive operation; durable, textured paint finish, in a wide selection of mar-resistant, easy-to-clean colors.

Any Base Cabinet can be made **MOBILE**.



TOP

C-1520-EP-1 FEATURES

One "B" power panel, one "B" drawer, one "C" drawer base cabinet.

Cabinet Dimensions: 15" wide, 33 3/4" high, 28" deep

"B" Power Panel supplied with EP-1 Electrical Kit

Dimensions: 15" wide, 6 7/32" high, 1" thick.

"B" Drawer

Drawer Face Dimensions: 15" wide, 6 11/16" high, 1" thick.

Drawer Body Inside Dimensions: 12" wide, 5" high, 25 1/2" deep.

Weight Capacity: 75 lbs.

"C" Drawer

Drawer Face Dimensions: 15" wide, 13 1/2" high, 1" thick.

Drawer Body Inside Dims.: 12" wide, 9 1/2" high, 25 1/2" deep.

Weight Capacity: 75 lbs.

All 15" wide "C" drawers accept standard hanging files.



FRONT



SIDE