


**BENCH HIGH BASE CABINET C-1542**

Teclab Base Cabinets are constructed of only top quality, heavy-duty, furniture-grade, cold-rolled steel, double-braced and welded at all critical joints; twin full-perimeter frames formed of heavy gauge steel, welded to provide strength and rigidity; heavy-duty drawers that glide effortlessly on precision steel ball bearings; cabinet doors that open and close with authority and precision using self-closing hinges; brushed aluminum pulls that are recessed for clean lines and unobtrusive operation; durable, textured paint finish, in a wide selection of mar-resistant, easy-to-clean colors.



**TOP**

**C-1542 FEATURES**

**One "A" drawer, one "E" door with LH pull and sliding bottom tray base cabinet.**

Cabinet Dimensions: 15" wide, 33 3/4" high, 28" deep

**"A" Drawer**

Drawer Face Dimensions: 15" wide, 3 1/4" high, 1" thick.

Drawer Body Inside Dimensions: 12" wide, 2" high, 25 1/2" deep.

Weight Capacity: 75 lbs.

*Top "A" drawer can be replaced with a writing surface.*


**"E" Door with Left Hand Pull**

Door Face Dimensions: 15" wide, 23 3/4" high, 1" thick.

**Sliding Tray**

Dimensions: 12" wide, 2" high, 25 1/2" deep.

Weight Capacity: 75 lbs.



**FRONT**



**SIDE**