


DESK HIGH SINK CABINET C-2005

Teclab Base Cabinets are constructed of only top quality, heavy-duty, furniture-grade, cold-rolled steel, double-braced and welded at all critical joints; twin full-perimeter frames formed of heavy gauge steel, welded to provide strength and rigidity; heavy-duty drawers that glide effortlessly on precision steel ball bearings; cabinet doors that open and close with authority and precision using self-closing hinges; brushed aluminum pulls that are recessed for clean lines and unobtrusive operation; durable, textured paint finish, in a wide selection of mar-resistant, easy-to-clean colors.



C-2005 FEATURES

One "B" panel, double "B" door base cabinet with sliding bottom tray.

Cabinet Dimensions: 30" wide, 27 ¾" high, 28" deep

"B" Panel

Dimensions: 30" wide, 6 11/16" high, 1" thick.

Top "A" drawer can be replaced with a writing surface.

"B" Doors

Door Face Dimensions: 14 3/8" wide, 16 15/16" high, 1" thick.

Sliding Tray

Dimensions: 27" wide, 2" high, 25 ½" deep.

Weight Capacity: 75 lbs.

