


## **PRODUCT DESCRIPTION**

# DESK HIGH BASE CABINET C-2103

Teclab Base Cabinets are constructed of only top quality, heavy-duty, furniture-grade, cold-rolled steel, double-braced and welded at all critical joints; twin full-perimeter frames formed of heavy gauge steel, welded to provide strength and rigidity; heavy-duty drawers that glide effortlessly on precision steel ball bearings; cabinet doors that open and close with authority and precision using self-closing hinges; brushed aluminum pulls that are recessed for clean lines and unobtrusive operation; durable, textured paint finish, in a wide selection of mar-resistant, easy-to-clean colors.





#### ΤΟΡ

### **C-2103 FEATURES**


**Double "E" door base cabinet with sliding bottom tray.** Cabinet Dimensions: 22" wide, 27 <sup>3</sup>/<sub>4</sub>" high, 28" deep.

"E" Doors

Door Face Dimensions: 11" wide, 23 <sup>3</sup>/<sub>4</sub>" high, 1" thick.

#### Sliding Tray

Dimensions: 19" wide, 2" high, 25 ½" deep. Weight Capacity: 75 lbs.



PD - C-2103 - v.1.3