


## STACKED BASE CABINET C-6013-30

Teclab Base Cabinets are constructed of only top quality, heavy-duty, furniture-grade, cold-rolled steel, double-braced and welded at all critical joints; twin full-perimeter frames formed of heavy gauge steel, welded to provide strength and rigidity; heavy-duty drawers that glide effortlessly on precision steel ball bearings; cabinet doors that open and close with authority and precision using self-closing hinges; brushed aluminum pulls that are recessed for clean lines and unobtrusive operation; durable, textured paint finish, in a wide selection of mar-resistant, easy-to-clean colors.

Any Base Cabinet can be made **MOBILE**.



### C-6013-30 FEATURES

**Four "A" drawers, two "B" drawers, two "C" drawers base cabinet.**  
Cabinet Dimensions: 30" wide, 61 1/2" high, 28" deep

#### "A" Drawer

Drawer Face Dimensions: 30" wide, 3 1/4" high, 1" thick.  
Drawer Body Inside Dimensions: 27" wide, 2" high, 25 1/2" deep.  
Weight Capacity: 75 lbs.

#### "B" Drawer

Drawer Face Dimensions: 30" wide, 6 11/16" high, 1" thick.  
Drawer Body Inside Dimensions: 27" wide, 5" high, 25 1/2" deep.  
Weight Capacity: 75 lbs.

#### "C" Drawer

Drawer Face Dimensions: 30" wide, 13 1/2" high, 1" thick.  
Drawer Body Inside Dimensions: 27" wide, 9 1/2" high, 25 1/2" deep.  
Weight Capacity: 75 lbs.

*All "C" drawers accept standard hanging files with cabinet accessory HFK-30.*

