

DESK HIGH BASE CABINET DL-1807

Teclab builds modular cabinetry of superior strength, thanks to twin full-perimeter frames. We also set door and drawer faces into the cabinet shell to protect the contents from spills.



DL-1807 FEATURES

Two "A" drawer, one "A" door with RH pull base cabinet.

Cabinet Dimensions: 18" wide, 28" high, 22" deep


"A" Drawer

Drawer Face Dimensions: 16 3/8" wide, 3 1/4" high, 1" thick.
 Drawer Body Inside Dims.: 13 7/8" wide, 2" high, 19 1/2" deep.
 Weight Capacity: 75 lbs.


Top "A" drawer can be replaced with a writing surface.

"A" Door with Right Hand Pull


Door Face Dimensions: 16 3/8" wide, 13 1/2" high, 1" thick.



TOP



FRONT



SIDE