

**DESK HIGH BASE CABINET DL-1808**

Teclab builds modular cabinetry of superior strength, thanks to twin full-perimeter frames. We also set door and drawer faces into the cabinet shell to protect the contents from spills.



**DL-1808 FEATURES**

**Two "A" drawer, one "A" door with LH pull base cabinet.**

Cabinet Dimensions: 18" wide, 28" high, 22" deep


**"A" Drawer**

Drawer Face Dimensions: 16 3/8" wide, 3 1/4" high, 1" thick.  
 Drawer Body Inside Dims.: 13 7/8" wide, 2" high, 19 1/2" deep.  
 Weight Capacity: 75 lbs.


*Top "A" drawer can be replaced with a writing surface.*

**"A" Door with Left Hand Pull**


Door Face Dimensions: 16 3/8" wide, 13 1/2" high, 1" thick.



**TOP**



**FRONT**



**SIDE**