

DESK HIGH BASE CABINET DL-1813

Teclab builds modular cabinetry of superior strength, thanks to twin full-perimeter frames. We also set door and drawer faces into the cabinet shell to protect the contents from spills.

Any Base Cabinet can be made **MOBILE**.

Height of a Desk High cabinet will increase by 3".



DL-1813 FEATURES

Five "A" drawer, one "B" drawer base cabinet.

Cabinet Dimensions: 18" wide, 28" high, 22" deep

"A" Drawer

Drawer Face Dimensions: 16 $\frac{3}{8}$ " wide, 3 $\frac{1}{4}$ " high, 1" thick.

Drawer Body Inside Dims.: 13 $\frac{7}{8}$ " wide, 2" high, 19 $\frac{1}{2}$ " deep.

Weight Capacity: 75 lbs.


Top "A" drawer can be replaced with a writing surface.

"B" Drawer


Drawer Face Dimensions: 16 $\frac{3}{8}$ " wide, 6 $\frac{1}{16}$ " high, 1" thick.

Drawer Body Inside Dims.: 13 $\frac{7}{8}$ " wide, 5" high, 19 $\frac{1}{2}$ " deep.


Weight Capacity: 75 lbs.



TOP



FRONT



SIDE