

**DESK HIGH BASE CABINET DL-1814**

Teclab builds modular cabinetry of superior strength, thanks to twin full-perimeter frames. We also set door and drawer faces into the cabinet shell to protect the contents from spills.



**DL-1814 FEATURES**

**Seven "A" drawer base cabinet.**

Cabinet Dimensions: 18" wide, 28" high, 22" deep

**"A" Drawer**

Drawer Face Dimensions: 16 3/8" wide, 3 1/4" high, 1" thick.  
 Drawer Body Inside Dimensions: 13 7/8" wide, 2" high, 20" deep.  
 Weight Capacity: 75 lbs.

*Top "A" drawer can be replaced with a writing surface.*

