

DESK HIGH BASE CABINET DL-2409

Teclab builds modular cabinetry of superior strength, thanks to twin full-perimeter frames. We also set door and drawer faces into the cabinet shell to protect the contents from spills.



DL-2409 FEATURES

Two "E" drawer base cabinet.

Cabinet Dimensions: 24" wide, 28" high, 22" deep

"E" Drawer

Drawer Face Dimensions: 22 3/8" wide, 11 13/16" high, 1" thick.
 Drawer Body Inside Dims.: 19 7/8" wide, 9 1/2" high, 19 1/2" deep.
 Weight Capacity: 75 lbs.

All "E" drawers accept standard hanging files with the addition of our Hanging File Adapter Kit.

