

DESK HIGH BASE CABINET DL-2410


Teclab builds modular cabinetry of superior strength, thanks to twin full-perimeter frames. We also set door and drawer faces into the cabinet shell to protect the contents from spills.



DL-2410 FEATURES

One "A" drawer, one "B" drawer, one "C" drawer base cabinet.

Cabinet Dimensions: 24" wide, 28" high, 22" deep



TOP

"A" Drawer

Drawer Face Dimensions: 22 3/8" wide, 3 1/4" high, 1" thick.
 Drawer Body Inside Dims.: 19 7/8" wide, 2" high, 19 1/2" deep.
 Weight Capacity: 75 lbs.

Top "A" drawer can be replaced with a writing surface.


"B" Drawer

Drawer Face Dimensions: 22 3/8" wide, 6 1/16" high, 1" thick.
 Drawer Body Inside Dims.: 19 7/8" wide, 5" high, 19 1/2" deep.
 Weight Capacity: 75 lbs.


"C" Drawer

Drawer Face Dimensions: 22 3/8" wide, 13 1/2" high, 1" thick.
 Drawer Body Inside Dims.: 19 7/8" wide, 9 1/2" high, 19 1/2" deep.
 Weight Capacity: 75 lbs.

All "C" drawers accept standard hanging files with the addition of our Hanging File Adapter Kit.



FRONT



SIDE