

**DESK HIGH BASE CABINET DL-2412**

Teclab builds modular cabinetry of superior strength, thanks to twin full-perimeter frames. We also set door and drawer faces into the cabinet shell to protect the contents from spills.



**DL-2412 FEATURES**

**Three "A" drawer, two "B" drawer base cabinet.**

Cabinet Dimensions: 24" wide, 28" high, 22" deep

**"A" Drawer**

Drawer Face Dimensions: 22 3/8" wide, 3 1/4" high, 1" thick.  
 Drawer Body Inside Dims.: 19 7/8" wide, 2" high, 19 1/2" deep.  
 Weight Capacity: 75 lbs.

*Top "A" drawer can be replaced with a writing surface.*

**"B" Drawer**

Drawer Face Dimensions: 22 3/8" wide, 6 1/16" high, 1" thick.  
 Drawer Body Inside Dims.: 19 7/8" wide, 5" high, 19 1/2" deep.  
 Weight Capacity: 50lbs.

