

DESK HIGH BASE CABINET DL-2416

Teclab builds modular cabinetry of superior strength, thanks to twin full-perimeter frames. We also set door and drawer faces into the cabinet shell to protect the contents from spills.

Any Base Cabinet can be made **MOBILE**.

Height of a Desk High cabinet will increase by 3".




DL-2416 FEATURES

One drawer base cabinet.


Cabinet Dimensions: 24" wide, 28" high, 22" deep

Single Drawer


Drawer Face Dimensions: 22 3/8" 23 3/4" high, 1" thick.
Weight Capacity: 75 lbs.



TOP



FRONT



SIDE

PD - DL-2416 - v.1.2