

DESK HIGH BASE CABINET DL-3003

Teclab builds modular cabinetry of superior strength, thanks to twin full-perimeter frames. We also set door and drawer faces into the cabinet shell to protect the contents from spills.



DL-3003 FEATURES

One "B" drawer, two "A" door base cabinet.

Cabinet Dimensions: 30" wide, 28" high, 22" deep

"B" Drawer

Drawer Face Dimensions: 28 3/8" wide, 6 1/16" high, 1" thick.
 Drawer Body Inside Dims.: 25 7/8" wide, 5" high, 19 1/2" deep.
 Weight Capacity: 75 lbs.

"A" Door

Door Face Dimensions: 14 1/4" wide, 13 1/2" high, 1" thick.

