

DESK HIGH BASE CABINET DL-3004

Teclab builds modular cabinetry of superior strength, thanks to twin full-perimeter frames. We also set door and drawer faces into the cabinet shell to protect the contents from spills.



DL-3004 FEATURES

Three "A" drawer, two "B" drawer base cabinet.

Cabinet Dimensions: 30" wide, 28" high, 22" deep

"A" Drawer

Drawer Face Dimensions: 28 3/8" wide, 3 1/4" high, 1" thick.
 Drawer Body Inside Dims.: 25 7/8" wide, 2" high, 19 1/2" deep.
 Weight Capacity: 75 lbs.

Top "A" drawer can be replaced with a writing surface.

"B" Drawer


Drawer Face Dimensions: 28 3/8" wide, 6 1/16" high, 1" thick.
 Drawer Body Inside Dims.: 25 7/8" wide, 5" high, 19 1/2" deep.
 Weight Capacity: 75 lbs.



TOP



FRONT



SIDE