

DESK HIGH BASE CABINET DL-3604

Teclab builds modular cabinetry of superior strength, thanks to twin full-perimeter frames. We also set door and drawer faces into the cabinet shell to protect the contents from spills.



DL-3604 FEATURES

Two half-width "A" drawer, six half-width "B" drawer base cabinet.

Cabinet Dimensions: 36" wide, 28" high, 22" deep

"A" Drawer (half-width)

Drawer Face Dimensions: 16 3/8" wide, 3 1/4" high, 1" thick.
 Drawer Body Inside Dims.: 13 7/8" wide, 2" high, 19 1/2" deep.
 Weight Capacity: 75 lbs.

Top "A" drawer can be replaced with a writing surface.

"B" Drawer (half-width)


Drawer Face Dims.: 16 3/8" wide, 6 1/16" high, 1" thick.
 Drawer Body Inside Dims.: 13 7/8" w, 5" h, 19 1/2" deep.
 Weight Capacity: 75 lbs.



TOP



FRONT



SIDE