

DESK HIGH SINK CABINET DL-3607

Teclab builds modular cabinetry of superior strength, thanks to twin full-perimeter frames. We also set door and drawer faces into the cabinet shell to protect the contents from spills.

DL-3607 FEATURES

One "B" panel, two "A" door base cabinet.

Cabinet Dimensions: 36" wide, 28" high, 22" deep

"B" Panel

Panel Face Dimensions: 34 $\frac{3}{8}$ " wide, 6 $\frac{1}{16}$ " high, 1" thick.

"A" Door

Door Face Dimensions: 17 $\frac{1}{4}$ " wide, 13 $\frac{1}{2}$ " high, 1" thick.

