

**BENCH HIGH BASE CABINET WL-1801**

Teclab builds modular cabinetry of superior strength, thanks to twin full-perimeter frames. We also set door and drawer faces into the cabinet shell to protect the contents from spills.




**WL-1801 FEATURES**

**One "D" door with RH pull base cabinet.**

Cabinet Dimensions: 18" wide, 35 1/2" high, 22" deep

**"D" Door with Right Hand Pull**


Door Face Dimensions: 16 3/8" wide, 27 3/16" high, 1" thick.



**TOP**



**FRONT**



**SIDE**

