

BENCH HIGH SINK CABINET WL-1808

Teclab builds modular cabinetry of superior strength, thanks to twin full-perimeter frames. We also set door and drawer faces into the cabinet shell to protect the contents from spills.

WL-1808 FEATURES

One "B" panel, One "C" door with LH pull base cabinet.


Cabinet Dimensions: 18" wide, 35 1/2" high, 22" deep

"B" Panel

Dimensions: 16 3/8" wide, 6 1/16" high, 1" thick.

"C" Door with Left Hand Pull


Door Face Dimensions: 16 3/8" wide, 20 3/8" high, 1" thick.



TOP



FRONT



SIDE