

## BENCH HIGH BASE CABINET WL-1815

Teclab builds modular cabinetry of superior strength, thanks to twin full-perimeter frames. We also set door and drawer faces into the cabinet shell to protect the contents from spills.



### WL-1815 FEATURES


#### Two "C" drawer base cabinet.

Cabinet Dimensions: 18" wide, 35 1/2" high, 22" deep

#### "C" Drawer

Drawer Face Dimensions: 16 3/8" wide, 13 1/2" high, 1" thick.  
Drawer Body Inside Dims.: 13 7/8" wide, 9 1/2" high, 19 1/2" deep.  
Weight Capacity: 75 lbs.


*All "C" drawers accept standard hanging files with the addition of our Hanging File Adapter Kit.*



TOP



FRONT



SIDE