

BENCH HIGH BASE CABINET WL-1824

Teclab builds modular cabinetry of superior strength, thanks to twin full-perimeter frames. We also set door and drawer faces into the cabinet shell to protect the contents from spills.




WL-1824 FEATURES

Open Base Cabinet with four sliding shelves.

Cabinet Dimensions: 18" wide, 35 1/2" high, 22" deep

Sliding Shelf


Dimensions: 13 7/8" wide, 2" high, 20" deep.
Weight Capacity: 75 lbs.



TOP



FRONT



SIDE

