

BENCH HIGH BASE CABINET WL-1826

Teclab builds modular cabinetry of superior strength, thanks to twin full-perimeter frames. We also set door and drawer faces into the cabinet shell to protect the contents from spills.



WL-1826 FEATURES

One Drawer Base Cabinet.


Cabinet Dimensions: 18" wide, 35 1/2" high, 22" deep

Single Drawer


Drawer Face Dimensions: 16 3/8" wide, 27 3/16" high, 1" thick.
Weight Capacity: 75 lbs.



TOP



FRONT



SIDE