

## BENCH HIGH BASE CABINET WL-2401

Teclab builds modular cabinetry of superior strength, thanks to twin full-perimeter frames. We also set door and drawer faces into the cabinet shell to protect the contents from spills.



### WL-2401 FEATURES

**One "D" door with RH pull base cabinet.**

Cabinet Dimensions: 24" wide, 35 1/2" high, 22" deep

**"D" Door with Right Hand Pull**

Door Face Dimensions: 22 3/8" wide, 27 3/16" high, 1" thick.

