

BENCH HIGH BASE CABINET WL-2402

Teclab builds modular cabinetry of superior strength, thanks to twin full-perimeter frames. We also set door and drawer faces into the cabinet shell to protect the contents from spills.




WL-2402 FEATURES

One "D" door with LH pull base cabinet.


Cabinet Dimensions: 24" wide, 35 1/2" high, 22" deep

"D" Door with Left Hand Pull


Door Face Dimensions: 22 3/8" wide, 27 3/16" high, 1" thick.



TOP



FRONT



SIDE