

BENCH HIGH BASE CABINET WL-2405

Teclab builds modular cabinetry of superior strength, thanks to twin full-perimeter frames. We also set door and drawer faces into the cabinet shell to protect the contents from spills.

Any Base Cabinet can be made **MOBILE**.



WL-2405 FEATURES

One "B" drawer, one "C" door with RH pull base cabinet.

Cabinet Dimensions: 24" wide, 35 1/2" high, 22" deep

"B" Drawer

Drawer Face Dimensions: 22 3/8" wide, 6 11/16" high, 1" thick.
Drawer Body Inside Dimensions: 19 7/8" w, 5" h, 19 1/2" deep.
Weight Capacity: 75 lbs.

"C" Door with Right Hand Pull

Door Face Dimensions: 22 3/8" wide, 20 3/8" high, 1" thick.

