

BENCH HIGH BASE CABINET WL-2426

Teclab builds modular cabinetry of superior strength, thanks to twin full-perimeter frames. We also set door and drawer faces into the cabinet shell to protect the contents from spills.



WL-2426 FEATURES

One drawer base cabinet.

Cabinet Dimensions: 24" wide, 35 1/2" high, 22" deep

Single Drawer

Drawer Face Dimensions: 22 3/8" wide, 27 3/16" high, 1" thick.
Weight Capacity: 75 lbs.

