

**BENCH HIGH BASE CABINET WL-3007**

Teclab builds modular cabinetry of superior strength, thanks to twin full-perimeter frames. We also set door and drawer faces into the cabinet shell to protect the contents from spills.



**WL-3007 FEATURES**

**Three "A" drawer, two "B" door base cabinet.**

Cabinet Dimensions: 30" wide, 35 1/2" high, 22" deep

**"A" Drawer**

Drawer Face Dimensions: 28 3/8" wide, 3 1/4" high, 1" thick.  
 Drawer Body Inside Dimensions: 25 7/8" w, 2" h, 19 1/2" deep.  
 Weight Capacity: 75 lbs.

**"B" Door**

Door Face Dimensions: 14 1/4" wide, 16 15/16" high, 1" thick.

