

**BENCH HIGH BASE CABINET WL-3610**

Teclab builds modular cabinetry of superior strength, thanks to twin full-perimeter frames. We also set door and drawer faces into the cabinet shell to protect the contents from spills.



**WL-3610 FEATURES**

**Four half width "A" drawer, two "C" door base cabinet.**


Cabinet Dimensions: 36" wide, 35 1/2" high, 22" deep

**"A" Drawer (half width)**


Drawer Face Dimensions: 16 3/8" wide, 3 1/4" high, 1" thick.  
 Drawer Body Inside Dims.: 13 7/8" w, 2" h, 19 1/2" deep.  
 Weight Capacity: 75 lbs.

**"C" Door**

Door Face Dimensions: 17 1/4" wide, 20 3/8" high, 1" thick.



**FRONT**



**SIDE**