

BENCH HIGH BASE CABINET WL-3613

Teclab builds modular cabinetry of superior strength, thanks to twin full-perimeter frames. We also set door and drawer faces into the cabinet shell to protect the contents from spills.

Any Base Cabinet can be made **MOBILE**.



WL-3613 FEATURES

Two half width "B" drawer, three full width "B" drawer base cabinet.


Cabinet Dimensions: 36" wide, 35 1/2" high, 22" deep

"B" Drawer (half width)

Drawer Face Dimensions: 16 1/4" w, 6 11/16" h, 1" thick.
Drawer Body Inside Dims.: 13 7/8" w, 5" h, 19 1/2" deep.
Weight Capacity: 75 lbs.

"B" Drawer (full width)


Drawer Face Dimensions: 34 3/8" wide, 6 11/16" high, 1" thick.
Drawer Body Inside Dims.: 31 7/8" w, 5" h, 19 1/2" deep.
Weight Capacity: 75 lbs.



TOP



FRONT



SIDE