

## BENCH HIGH BASE CABINET WL-3625

Teclab builds modular cabinetry of superior strength, thanks to twin full-perimeter frames. We also set door and drawer faces into the cabinet shell to protect the contents from spills.

Any Base Cabinet can be made **MOBILE**.



### WL-3625 FEATURES

#### Two full width "C" drawer base cabinet.

Cabinet Dimensions: 36" wide, 35 1/2" high, 22" deep

#### "C" Drawer (full width)

Drawer Face Dimensions: 34 3/8" w, 13 1/2" h, 1" thick.

Drawer Body Inside Dims.: 31 7/8" w, 9 1/2" h, 19 1/2" deep.

Weight Capacity: 75 lbs.

*All "C" drawers accept standard hanging files with the addition of our Hanging File Adapter Kit.*

