

BENCH HIGH BASE CABINET WL-3653

Teclab builds modular cabinetry of superior strength, thanks to twin full-perimeter frames. We also set door and drawer faces into the cabinet shell to protect the contents from spills.

Any Base Cabinet can be made **MOBILE**.




WL-3653 FEATURES

Open base cabinet with four sliding shelves.


Cabinet Dimensions: 36" wide, 35 1/2" high, 22" deep

Sliding Shelf

Dimensions: 31 7/8" wide, 2" high, 19 1/2" deep.
Weight Capacity: 75 lbs.



TOP



FRONT



SIDE