

BENCH HIGH BASE CABINET WL-3654

Teclab builds modular cabinetry of superior strength, thanks to twin full-perimeter frames. We also set door and drawer faces into the cabinet shell to protect the contents from spills.

Any Base Cabinet can be made **MOBILE**.



WL-3654 FEATURES

Two drawer base cabinet.

Cabinet Dimensions: 36" wide, 35 1/2" high, 22" deep

Full Height Drawer

Drawer Face Dimensions: 16 3/8" wide, 27 3/16" high, 1" thick.



TOP



FRONT



SIDE