

BENCH HIGH BASE CABINET WL-4249

Teclab builds modular cabinetry of superior strength, thanks to twin full-perimeter frames. We also set door and drawer faces into the cabinet shell to protect the contents from spills.



WL-4249 FEATURES


Open base cabinet with four sliding shelves.

Cabinet Dimensions: 42" wide, 35 1/2" high, 22" deep


Sliding Shelf

Dimensions: 37 7/8" wide, 2" high, 19 1/2" deep.

Weight Capacity: 75 lbs.



TOP



FRONT



SIDE