

## ISLAND BASE CABINET WL-4604-1

Teclab builds modular cabinetry of superior strength, thanks to twin full-perimeter frames. We also set door and drawer faces into the cabinet shell to protect the contents from spills.

### WL-4604 FEATURES

#### Four "B" drawer base cabinet.

Cabinet Dimensions: 46" wide, 35 1/2" high, 22" deep

#### "B" Drawer

Drawer Face Dimensions: 44 3/8" wide, 6 1/16" high, 1" thick.

Dwr. Body Inside Dims.: 41 7/8" wide, 5" high, 19 1/2" deep.

Weight Capacity: 75 lbs.

#### NOTE:

Island Cabinets are also available with standard full-width bases for use when not positioned at the end of an island.

