

**BENCH HIGH BASE CABINET WL-4818**

Teclab builds modular cabinetry of superior strength, thanks to twin full-perimeter frames. We also set door and drawer faces into the cabinet shell to protect the contents from spills.



**WL-4818 FEATURES**

**Two "A" drawer, three "B" drawer base cabinet.**

Cabinet Dimensions: 48" wide, 35 1/2" high, 22" deep

**"A" Drawer (full width)**

Drawer Face Dimensions: 46 3/8" wide, 3 1/4" high, 1" thick.  
 Drawer Body Inside Dims.: 43 7/8" w, 2" h, 19 1/2" deep.  
 Weight Capacity: 75 lbs.


*Top "A" drawer can be replaced with a writing surface.*

**"B" Drawer (full width)**

Drawer Face Dimensions: 46 3/8" w, 6 1/16" h, 1" thick.  
 Drawer Body Inside Dims.: 43 7/8" w, 5" h, 19 1/2" deep.  
 Weight Capacity: 75 lbs.



**TOP**



**FRONT**



**SIDE**