

## ISLAND BASE CABINET WL-5802-1

Teclab builds modular cabinetry of superior strength, thanks to twin full-perimeter frames. We also set door and drawer faces into the cabinet shell to protect the contents from spills.

### WL-5802 FEATURES

**Two "B" drawer, two "C" door base cabinet.**

Cabinet Dimensions: 58" wide, 35 1/2" high, 22" deep

**NOTE:**

Island Cabinets are also available with standard full-width bases for use when not positioned at the end of an island.

**"B" Drawer**

Drawer Face Dimensions: 22 3/8" wide, 6 11/16" high, 1" thick.  
 Drawer Body Inside Dims.: 19 7/8" wide, 5" high, 19 1/2" deep.  
 Weight Capacity: 75 lbs.

**Center Panel**

Dimensions: 11 3/8" wide, 27 3/16" high, 1" thick.

**"C" Door**

Door Face Dimensions: 22 3/8" wide, 20 3/8" high, 1" thick.

