


## DRAWER DIVIDER KIT DKA-2201

Organize your Teclab cabinets with these drop-in kits. Made to fit “A” and “B” drawers in all cabinet widths. Each kit comes with 2 side rails, a drawer mat, and enough spacers to configure drawers in many ways.

### DKA-2201 Divider Kit for 22” wide “A” drawer

Configuration: 10 spaces (9”wide)



**Dividers:** (stainless steel)  
8 (9.47” x 1.75”)

**Rails:** (finished in same manner as drawer)  
1 Left Side  
1 Right Side  
1 Center

**Drawer Mats:**  
2 (9.47”wide)