


DRAWER DIVIDER KIT DKA-3006

Organize your Teclab cabinets with these drop-in kits. Made to fit “A” and “B” drawers in all cabinet widths. Each kit comes with 2 side rails, a drawer mat, and enough spacers to configure drawers in many ways.

DKA-3006 Divider Kit for 30” wide “A” drawer

Configuration: 18 spaces (7.8” wide)



Dividers: (stainless steel)
15 (8.223” x 1.75”)

Rails: (finished in same manner as drawer)
1 Left Side
1 Right Side
2 Centers

Drawer Mats:
3 (8.223” wide)