



## DRAWER DIVIDER KIT DKB-3608

Organize your Teclab cabinets with these drop-in kits. Made to fit “A” and “B” drawers in all cabinet widths. Each kit comes with 2 side rails, a drawer mat, and enough spacers to configure drawers in many ways.

### DKB-3608 Divider Kit for 36” wide “B” drawer

Configuration: 40 spaces (7.2” wide)



**Dividers:** (stainless steel)  
36 (7.638” x 3.5”)

**Rails:** (finished in same manner as drawer)  
1 Left Side  
1 Right Side  
3 Centers

**Drawer Mats:**  
4 (7.638” wide)