


## DRAWER DIVIDER KIT WDKA-1801

Organize your Teclab cabinets with these drop-in kits. Made to fit “A” and “B” drawers in all cabinet widths. Each kit comes with 2 side rails, a drawer mat, and enough spacers to configure drawers in many ways.

### WDKA-1801 Divider Kit for 18” wide “A” drawer

Configuration: 4 spaces (13.1” wide)



**Dividers:** (Stainless Steel)  
3 (13.556” x 1.75”)

**Rails:** (finished in same manner as drawer)  
1 Left Side  
1 Right Side  
1 Center

**Drawer Mat:**  
1 (13.556” wide)