


DRAWER DIVIDER KIT WDKA-2103

Organize your Teclab cabinets with these drop-in kits. Made to fit “A” and “B” drawers in all cabinet widths. Each kit comes with 2 side rails, a drawer mat, and enough spacers to configure drawers in many ways.

WDKA-2103 Divider Kit for 21” wide “A” drawer

Configuration: 10 spaces (7.7” wide)



Dividers: (Stainless Steel)
8 (8.141” x 1.75”)

Rails: (finished in same manner as drawer)
1 Left Side
1 Right Side
1 Center

Drawer Mats:
2 (8.141” wide)