



DRAWER DIVIDER KIT WDKA-3604

Organize your Teclab cabinets with these drop-in kits. Made to fit "A" and "B" drawers in all cabinet widths. Each kit comes with 2 side rails, a drawer mat, and enough spacers to configure drawers in many ways.

WDKA-3604 Divider Kit for 36" wide "A" drawer

Configuration: 16 spaces (7.4" wide)



Dividers: (Stainless Steel)
12 (7.802" x 1.75")

Rails: (finished in same manner as drawer)
1 Left Side
1 Right Side
3 Centers

Drawer Mats:
4 (7.802" wide)