

LABORATORY STOOL

Choose our Laboratory Stool for top performance in a wide variety of environments including: industry, education, laboratory, and health care. This stool comes standard with a 5-legged tubular base/foot ring and single-wheel ball bearing casters for hard floors.

STANDARD FEATURES

Base:

Rugged 5-legged tubular steel base with 20" diameter attached footring, 23" leg spread.

Seat:

Round, 16" diameter, upholstered, with 4" of high-resilience molded foam cushioning.

Seat Height:

CLS-1726 is adjustable from 17" to 26".
CLS-2434 is adjustable from 24" to 34".

OPTIONS

Glides:

Nylon glides.

Metal Parts:

A black powder-coated finish comes standard.

A chrome-plated finish is available as an option.

Upholstery:

Available in gray or black. Other colors available on request.



CLS-1726
(shown)

CLS-2434