

Teclab has a select line of laboratory seating built to the highest standards in the industry. All seating is ergonomically designed for outstanding comfort, total body support and correct seating posture. All things leading to increased worker productivity. Each chair can be configured to your unique laboratory environment allowing you to specify such items as seating materials, chair heights and arm rests. Virtually anything to make the most comfortable seating possible.

Contact us and we'll discuss all the options available. We can help you make your next chair selection a vital addition to your laboratory while also making the most economic sense.



CHR-2432-ESD
(shown)

CHR-2432-SP
(Cleanroom)

PREMIUM HIGH RING ESD CHAIR & CLASS 100 CLEANROOM CHAIRS

These chairs come with a synchronized Fluid-Motion mechanism, which allows the seat and back to move in unison, providing continuous lumbar support to the user. Backrest movement "free floats" or can be locked at any point in the range. A sensitive tension control allows easy adjustment of the backrest pressure. Integral lumbar support adjusts to ensure correct seated posture. Comes standard with a 20" height adjustable foot ring, composite 5-star base and rugged 2" dual-wheel ESD casters.



CDC-1622

DESK HIGH CHAIR

This chair comes standard with a seat and back that adjust independently. Backrest angle follows the user's movement by pneumatic cylinder, and can be locked at any point in the range. Integral lumbar support adjusts to ensure correct posture. Comes standard with a composite 5-star base and rugged 2" dual-wheel casters.



CHR-2432

PREMIUM QUALITY HIGH RING CHAIR

This chair comes standard with a seat and back that adjust independently. Backrest angle follows the user's movement by pneumatic cylinder, and can be locked at any point in the range. Integral lumbar support adjusts to ensure correct posture. Comes standard with a 20" height adjustable foot ring, composite 5-star base and rugged 2" dual-wheel casters.



CHR1-2432

SEATWISE HIGH RING CHAIR

Seatwise intermediate and tall height chairs are equipped with a wide backrest that supports the user's back as they move. This feature is equipped with a tension adjustment and variable back lock. Fingertip adjustment of the seat height and backrest height adjustment with integral lumbar support make this an extremely comfortable, affordable ergonomic chair. Comes standard with a 20" diameter height adjustable foot ring, composite 5-star base and rugged 2" dual-wheel casters.



CHM-1822

MODERN MESH DESK CHAIR

A contemporary mesh chair designed to create the ultimate experience in comfort, support, and durability. The modern design provides features and functionality for today's work environment.

Mesh Chair, Contoured Backrest With Lumbar Support, Adjustable Arms, Black Nylon Base, All-Purpose Casters, 18.75" - 22.25" Seat Height Adj., 275lb Weight Capacity, 5yr Warranty.



CLS-1726
(desk-height)

CLS-2432
(bench-height)

LABORATORY STOOL

Choose our Laboratory Stool for top performance in a wide variety of environments including: industry, education, laboratory, and health care. This stool comes standard with a polished aluminum base and casters for hard floors.



CHU-2432

PREMIUM QUALITY HIGH RING POLYURETHANE CHAIR

This high ring chair has durable, contoured polyurethane seat and back with ventilated ribs for lasting comfort and support. Pneumatic seat height adjustment, backrest angle and height adjustment are standard. Chair comes equipped with dual wheel rubber casters for hard floors or plastic casters for carpeted areas. Sturdy five-star reinforced plastic base with an 18" diameter adjustable chrome foot ring also standard.