

**DECK MOUNTED GAS, AIR OR VACUUM
DOUBLE NOZZLE (GAV-D-DN)**

FEATURES:

- Straight Serrated Nozzle
- Metal Lever Handle with air, gas and vac service buttons
- Service: Natural Gas, Air or Vacuum
- Turret with Two Side Outlets @ 180° and Shank
- Solid Brass Valve Construction
- Fully Assembled and Tested

PERFORMANCE SPECIFICATIONS

- Rated Operating Pressure: 20-125 PSI
- Rated Operating Temperature: 40-140° F

COMPLIANT TO:

- ANSI Z21.5a
- CGA 9.1M

