

**WALL MOUNTED GAS, AIR OR VACUUM  
SINGLE NOZZLE (GAV-W-SN)**

**FEATURES:**

- Straight Serrated Nozzle
- Metal Lever Handle with air, gas and vac service buttons
- Service: Natural Gas, Air or Vacuum
- Straight Wall Flange With Shank
- Solid Brass Valve Construction
- Fully Assembled and Tested

**PERFORMANCE SPECIFICATIONS**

- Rated Operating Pressure: 20-125 PSI
- Rated Operating Temperature: 40-140° F

**COMPLIANT TO:**

- ANSI Z21.5a
- CGA 9.1M

